

Supplementary Information

Figure S1. Control experiments of immunostaining of *S. domuncula* tissue cross-sections. Sections were treated with blocking solution (PBS/BSA 0.3%, 1 h RT) and Cy3-labelled secondary antibody (1 h, RT). (a) Light microscope image and (b) corresponding fluorescent image showing no signal. Scale bar: 100 μm .




Figure S2. Light microscope images of *S. domuncula* tissue cross-sections after staining with Perls' Prussian Blue solution. Small spheres dispersed throughout tissue cross-sections are clearly observed.

