

Electrostatic surface properties of blood and semen extracellular vesicles: Implications of sialylation and HIV-induced changes on EV internalization

Hussein Kaddour¹, Tyler D. Panzner¹, Jennifer L. Welch^{2,3,4}, Nadia Shouman¹, Mahesh Mohan⁵, Jack T. Stapleton^{2,3,4}, and Chioma M. Okeoma^{1,*}

¹ Department of Pharmacology, Stony Brook University Renaissance School of Medicine, Stony Brook, NY 11794-8651, USA

² Department of Microbiology and Immunology, Carver College of Medicine, University of Iowa, Iowa City, Iowa, USA

³ Medical Service, Iowa City Veterans Affairs Medical Center, Iowa City, IA 52246-2209, USA

⁴ Department of Internal Medicine, Carver College of Medicine, University of Iowa, Iowa City, IA 52242-1109, USA

⁵ Host Pathogen Interaction Program, Southwest National Primate Research Center, Texas Biomedical Research Institute, San Antonio, Texas, 78227, USA; mmohan@txbiomed.org

* Correspondence: chioma.okeoma@stonybrook.edu; Tel.: +1-631-444-3014

Supplementary file

This file contains:

- Table S1. Parameters of the ζ -potential measurements with ZetaView.
- Figure S1. Particle size measurements of HIV- and HIV+ BEVs and SEVs.
- Figure S2. Representative plots of ζ -potential measurements for HIV- BEVs and SEVs.
- Figure S3. Representative plots of ζ -potential measurements for HIV+ BEVs and SEVs.

Table S1. Parameters of the ζ -potential measurements with Zeta-View.

Donor ID	EV	Dilution factor	particles per frame [¥]	number of traced particles [£]
1	BEV	160,000	151 ± 5	2128 ± 56
	SEV	160,000	123 ± 3	2895 ± 65
2	BEV	200,000	50 ± 14	1337 ± 335
	SEV	200,000	44 ± 5	753 ± 247
3	BEV	160,000	49 ± 9	1687 ± 198
	SEV	40,000	56 ± 14	1943 ± 380
4	BEV	80,000	73 ± 19	1492 ± 202
	SEV	80,000	49 ± 24	1559 ± 173
5	BEV	40,000	180 ± 97	1824 ± 384
	SEV	80,000	149 ± 66	3576 ± 620
6	BEV	80,000	169 ± 8	4241 ± 215
	SEV	80,000	97 ± 22	3129 ± 393
7	BEV	160,000	66 ± 49	2015 ± 284
	SEV	160,000	166 ± 65	3487 ± 598
8	BEV	160,000	100 ± 31	3489 ± 423
	SEV	80,000	127 ± 80	2668 ± 533
9	BEV	20,000	149 ± 9	3634 ± 72
	SEV	160,000	93 ±	3117 ± 304
10	BEV	40,000	146 ±	3829 ± 360
	SEV	80,000	106 ±	2118 ± 318
11	BEV	80,000	49 ± 37	2546 ± 642
	SEV	40,000	120 ± 14	3118 ± 452
12	BEV	200,000	175 ± 9	4602 ± 107
	SEV	40,000	94 ± 24	3279 ± 276
13	BEV	80,000	45 ± 6	1730 ± 318
	SEV	40,000	130 ± 11	3695 ± 259
1+	BEV	40,000	241 ± 30	4874 ± 470
	SEV	80,000	169 ± 7	4490 ± 100
2+	BEV	40,000	117 ± 23	3910 ± 334
	SEV	40,000	241 ± 18	4654 ± 420
3+	BEV	40,000	98 ± 10	3231 ± 240
	SEV	40,000	98 ± 19	3430 ± 251
4+	BEV	160,000	66 ± 10	1382 ± 240
	SEV	40,000	41 ± 2	1378 ± 115
5+	BEV	20,000	112 ± 27	2761 ± 360
	SEV	20,000	48 ± 1	1305 ± 119
6+	BEV	20,000	62 ± 17	2057 ± 216
	SEV	20,000	41 ± 1	1093 ± 86
7+	BEV	40,000	70 ± 9	2661 ± 834
	SEV	80,000	70 ± 32	1986 ± 413
8+	BEV	40,000	49 ± 3	1757 ± 54
	SEV	40,000	80 ± 2	2389 ± 3
9+	BEV	40,000	80 ± 13	2007 ± 236
	SEV	40,000	125 ± 10	1021 ± 105
10+	BEV	80,000	80 ± 7	2883 ± 85
	SEV	160,000	75 ± 12	2684 ± 138
11+	BEV	80,000	43 ± 25	1459 ± 960
	SEV	80,000	55 ± 4	2177 ± 164
12+	BEV	80,000	72 ± 28	1796 ± 94
	SEV	80,000	39 ± 6	1409 ± 845
13+	BEV	160,000	89 ± 10	2714 ± 74
	SEV	80,000	42 ± 19	1752 ± 100

¥ this parameter shows if samples were diluted to the right range where NTA is most accurate. Good PPF range is 40 – 250

£ this parameter shows the total number of particles analyzed during a measurement

* Errors represent standard deviation of 3 to 10 measurements per sample. Higher number of measurements was taken to test the stability of certain samples, which turned to be stable

Bold = samples from matched donors




Figure S1. Particle size measurements of HIV- and HIV+ BEVs and SEVs. Purified EVs were diluted to the appropriate concentration with ultrapure H₂O (1:320,000 – 1:20,000) and size measurements (1 cycle at 11 positions) were undertaken using ZetaView PMX 110 (Particle Metrix, Mebane, NC, USA) and the corresponding software (ZetaView v8.5.5.2). Post-acquisition parameters were fixed as follows: min brightness, 20; min size, 10; max size, 4000; and tracelength, 15. Camera control settings were as follows: sensitivity, 90; frame rate, 30; and shutter, 70. (a) HIV- BEVs and SEVs. n=13, 10 unmatched- and 3 matched-donor samples. The number above the bars represents the *p*-value of an unpaired t-test with Welch's correction. (b) HIV+ BEVs and SEVs. n=13 matched-donor samples. The number above the bars represents the *p*-value of a paired and parametric t-test. (c,d) Comparison of HIV- and HIV+ BEVs (c) and SEVs (d). Numbers above bars represent *p*-values of unpaired t-tests with Welch's correction. For all graphs, each data point represents mean size of triplicate measurements, numbers inside bars correspond to the mean of the means, and error bars correspond to S.D.




Figure S2. Representative plots of z-potential measurements for HIV- samples. (a) BEVs and (b) SEVs, from 10 non-autologous donors (D1-D10). (c) BEVs and SEVs from 3 matched donors (D11-13). Red dashed vertical line delimits 0 mV. Numbers inside graphs are mean values \pm S.D. of quintuplet measurements.




Figure S3. Representative plots of z-potential measurements for HIV+ samples. (a) BEVs and (b) SEVs from 13 autologous donors (D1-D13). Red dashed vertical line delimits 0 mV. Numbers inside graphs are mean values \pm S.D. of quintuplet measurements.