



Nanobody Technology for Mycotoxin Detection: Current Status and Prospects

Ting He¹, Jiang Zhu¹, Yao Nie¹, Rui Hu¹, Ting Wang², Peiwu Li^{2,*}, Qi Zhang^{2,*} and Yunhuang Yang^{1,*}

- State Key Laboratory of Magnetic Resonance and Atomic Molecular Physics, National Center for Magnetic Resonance in Wuhan, Wuhan Institute of Physics and Mathematics, Chinese Academy of Sciences, Wuhan 430071, China; ht210@wipm.ac.cn (T.H.); jiangzhu@wipm.ac.cn (J.Z.); nieyao@wipm.ac.cn (Y.N.); hurui@wipm.ac.cn (R.H.)
- ² Oil Crops Research Institute of the Chinese Academy of Agricultural Sciences, Wuhan 430062, China; wangting963@163.com
- * Correspondence: peiwuli@oilcrops.cn (P.L.); zhangqi01@caas.cn (Q.Z.); yang_yh@wipm.ac.cn (Y.Y.); Tel.: +86-27-8681-2943 (P.L. & Q.Z.); +86-27-8719-9541 (Y.Y.)

Received: 10 April 2018; Accepted: 27 April 2018; Published: 29 April 2018



Abstract: Mycotoxins, which are toxic, carcinogenic, and/or teratogenic, have posed a threat to food safety and public health. Sensitive and effective determination technologies for mycotoxin surveillance are required. Immunoassays have been regarded as useful supplements to chromatographic techniques. However, conventional antibodies involved in immunoassays are difficult to be expressed recombinantly and are susceptible to harsh environments. Nanobodies (or VHH antibodies) are antigen-binding sites of the heavy-chain antibodies produced from Camelidae. They are found to be expressed easily in prokaryotic or eukaryotic expression systems, more robust in extreme conditions, and facile to be used as surrogates for artificial antigens. These properties make them the promising and environmentally friendly immunoreagents in the next generation of immunoassays. This review briefly describes the latest developments in the area of nanobodies used in mycotoxin detection. Moreover, by integrating the introduction of the principle of nanobodies production and the critical assessment of their performance, this paper also proposes the prospect of nanobodies in the field of food safety in the foreseeable future.

Keywords: Nanobody; Mycotoxins; Immunoassay; Food safety

Key Contribution: This review thoroughly summarizes the advances and progress of nanobody technology applied to detection of mycotoxin contamination in agricultural products. It reviews the critical control points of nanobody production; evaluates their performance in mycotoxin detection; and discusses the prospects of nanobody technology in the future.

1. Introduction

Mycotoxins are toxic metabolites produced by fungal species that have adverse effects on humans, animals, and crops, resulting in significant medical costs caused by mycotoxin-based diseases and large economic losses in international trade [1]. Aflatoxins (AFs), ochratoxins (OTs), zearelenone (ZEN), fumonisins (FBs), deoxynivalenol (DON), and other trichothecenes are the major mycotoxins of greatest agro-economic importance [2]. In light of the wide spread of mycotoxins, many countries set action regulations for mycotoxins in various food and foodstuff. Increasingly sensitive and effective determination technologies for mycotoxins surveillance are required. Table 1 displays some predominant mycotoxins and their major fungal producers, contaminated hosts, toxic effects [3], and their maximum permitted levels in foods according to European legislation.



There have been many efforts related to mycotoxin determination, typically using highperformance liquid chromatography (HPLC) [4], liquid chromatography-tandem mass spectrometry (LC-MS/MS) [5], and immunoassays [6,7]. Among these methods, immunoassays have been regarded as useful supplements to chromatographic techniques due to their merits in sensitivity, specificity, and low cost [8]. Large molecules can be easily determined by a sandwich-format immunoassay [9,10]. Small analytes, such as mycotoxins, are not large enough to be simultaneously recognized by two antibodies, so the competitive immunoassay formats are normally used, which need a hapten to compete with the target analyte for binding to the antibody [11].

Mycotoxins	Major Producers	Major Host	Toxin Effects	Maximum Level
Aflatoxins (AFB ₁ , AFB ₂ , AFG ₁ , AFG ₂)	Aspergillus flavus, A. parasiticus	maize, nuts, cottonseed, peanuts	hepatotoxicity, cancer, immunosuppression	8.0 μg/kg (AFB ₁); 15 μg/kg (total aflatoxins) (Commission Regulation, EU No. 165/2010) [12]
Fumonisins (FB ₁ ,FB ₂)	Fusarium verticillioides, F. proliferatum	maize	hepatotoxicity, cancer, pulmonary edema, leukoencephalomalacia	4000 μg/kg in maize; 1000 μg/kg in maize-based foods (EC No. 1126/2007) [13]
Deoxynivalenol (DON)	F.graminearum, F. culmorum	maize, wheat, barley	gastrointestinal toxicity, inflammation of central nervous system	1250 μg/kg in cereals (EC No. 1126/2007)
Zearalenone (ZEN)	F. culmorum, F.graminearum	maize, barley, wheat, rice	infertility, abortion	350 μg/kg in maize; 100 μg/kg in other cereals (EC No. 1126/2007)
OchratoxinA (OTA)	Penicillium verrucosum, A.ochraceus	cereal-derived products, spices, wine	nephrotoxic, carcinogenic, teratogenic, immunotoxic effects	5.0 μg/kg in maize (EC No. 1881/2006)
Citrinin (CIT)	P.citrinum, P. verruscosum, P.expansum, A.terreus, Monascus ruber	corn, wheat, barley, rice	nephrotoxic, hepatotoxic, immunotoxic, carcinogenic effects	2000 μg/kg (EU No. 212/2014) [14]

Table 1. The major fungal producers, hosts, toxic effects of mycotoxins, and their maximum permitted levels in food according to European legislation.

A large number of antibodies against mycotoxins have been developed, including polyclonal antibodies (pAbs) [15,16], monoclonal antibodies (mAbs) [17,18], and recombinant antibodies (rAbs) [19,20]. Since heavy-chain-only antibodies (HcAbs) that are naturally devoid of light chain and lacking the CH1 domains were first discovered in the serum of the camel (*Camelus dromedarius*) [21], camelid antibodies have gained increasing attention due to their convenience in biotechnological applications. While conventional mAbs and pAbs are composed of two heavy and two light chains with variable domains of heavy chain (VH) and light chain (VL) contributing to form a paratope of the antibody, the camelid HcAbs are composed of only heavy chains, of which the antigen-binding sites are formed only by the variable domain of the heavy chain (VHH) [22]. These VHHs can be expressed recombinantly and have been termed VHH antibodies, or single domain antibodies (sdAbs) [23]. Due to the nanometer-scale size ($4 \text{ nm} \times 2.5 \text{ nm} \times 3 \text{ nm}$) of the VHH antibodies, they are also referred to as nanobodies (Nbs) by a biopharmaceutical company, Belgium Ablynx [24]. The representative structures of conventional antibody and nanobody derived from camelid HcAb are diagramed in Figure 1.

Sequence analysis and elucidation of the crystal structure have revealed several structural features of Nb fragments. The single domain nature of Nbs provides many advantages over conventional antibodies and other recombinant antibody fragments such as antigen-binding fragment (Fab), variable fragment (Fv), and single-chain variable fragment (scFv) [25–27]. The most characteristic feature of Nbs is the presence of amino acid substitutions in the framework-two region (FR2), making Nbs more soluble than conventional antibody fragments [28]. The extended flexible complementarity-

determining region 3 (CDR3) makes Nbs capable of recognizing antigenic sites that are normally inaccessible or cryptic for pAbs and mAbs such as enzyme active sites and cryptic viral epitopes [29,30]. In particular, their small size (12~15 kDa), ease of cloning and expressing in various expression systems, high solubility and stability, resistance to proteolysis, and facile genetic manipulation, make Nbs an excellent alternative to producing immunoreagents for immunoassays [31]. For example, Nbs can be used in biosensors with higher immobilizing densities [32], used for affinity purification medication [33], used in the design of anti-idiotypic antibodies [34] and as modular building units for multivalent or multifunctional constructs [35].



Figure 1. Schematic representation of a (A) conventional antibody (IgG) and (B) heavy chain camelid antibody. The isolated variable domain of the latter is called nanobody.

Nanobody technology has been initially applied in diagnostic and therapeutic areas [36–38]. To date, most reported biotechnological applications of Nbs refer to proteins or other macromolecular antigens, while rare reports are related to Nbs used against small hapten (molecular weight below 1500 Da) for the reason of the sensitivities of most assays at the μ M range that are somewhat unfavorable for practical application [39–42]. However, a growing number of studies have demonstrated that highly-sensitive Nbs maybe successfully isolated by sophisticated design of the selection procedure from the small subpopulation of high-affinity Nbs pool. There have been some successful examples regarding the use of Nbs to effectively recognize small molecules. Those compounds, such as the mycotoxins AFB₁, OTA, and 15-acetyl-DON, and some environmental contaminants, were reviewed by Bever et al. [43].

Recently, there have been various forms of Nbs used in mycotoxin immunodetections. In addition to the hapten-specific Nbs commonly used as primary antibodies involved in immunoassay, anti-idiotypic Nbs (AI-Nbs), which are generated from mAbs-immunized camelid and which are capable of targeting the antigenic determinants of the primary antibodies, can be used for immunoassays as surrogates for competing antigens and can significantly improve the sensitivity of the assays [44]. For many analytical and diagnostic applications, Nbs have been directly used as phage-borne antibody fragments, relying on the multiple copies of M13 phage particles to increase the assay sensitivity [45]. Likewise, the phage DNA can be used for ultrasensitive detection of the phage-borne antibody by utilizing a combination of real-time polymerase chain reaction (RT-PCR) in different immunoassay formats, which widens the dynamic range of the assay by several orders of magnitude [46,47].

Considering the predominant traits of Nbs, they could meet the field-testing requirement of reagents with greatly enhanced stability for analytical methods and become promising alternative immunoreagents for mycotoxin detection. However, few reviews have provided detailed summaries of the development of the Nb-based immunoassays utilized for mycotoxin determination. This review focuses on the introduction of Nb technology for mycotoxin detection in the field of food safety and refers to recently reported studies, integrating the introduction of the principle of Nb development

and providing a detailed description of their performance. This paper also proposes the prospect of Nbs in the area of analytical chemistry in the foreseeable future.

2. Principle of Developing Chemical-Specific Nb

Recombinant antibody techniques were used to create phage-displayed Nb libraries. To take advantage of the exquisite in vivo antibody-affinity maturation process, animals are often immunized prior to antibody isolation. It was noted that immune libraries lead more directly to anti-hapten Nbs with higher affinity and sensitivity than that of nonimmune libraries [43], although a few Nbs against small molecules (e.g., picloram [39] and auxin [40]) and AI-Nbs have been selected from native library [48–50].

Members of the Camelid family, including dromedary, llamas, and alpaca, are easily immunized following standard procedures [51]. Small haptens are compounds with low molecular weight and low immunogenic, customarily conjugated to a carrier protein (e.g., bovine serum albumin, BSA; ovalbumin, OVA; keyhole limpet hemocyanin, KLH, etc.) to induce an immune response. After immunization with weekly boosts, the total RNA is extracted from peripheral blood lymphocytes isolated from whole blood, and its cDNA synthesized using reverse transcriptase techniques [52]. The Nb-encoding fragments from IgG2 and IgG3-a subset of IgG antibodies lacking light chains are amplified and cloned into M13 phagemid vector to construct the phage-displayed Nb library. The resulting phage-Nb library can be used for panning target-specific binders.

The panning strategy of improving selection pressure has been proved to be valid in isolating positive phage-Nb clones with favorable affinity and sensitivity [53–55]. Non-specific phage-Nbs were minimized by increasingly stringent washings and desorption of phage input/output to the carrier protein. Moreover, blocking reagents should be switched alternatively in each panning process which is beneficial to minimize the nonspecific absorption. There are two regular methods for the elution of phage binding to the immobilized target in the process of biopanning. One is polar elution (for example, 0.2 M glycine-HCl pH 2.2; or 0.1 M triethylamine pH 11.0 for elution of phage-Nb), which is frequently used for a isolation of Nbs against large molecular antigens or in the case of lack free competitors [48,56]. The other is competitive elution, which is mostly applied to a selection of Nbs for small molecules, such as mycotoxins and chemical contaminants [53,55,57]. The construction of phage-displayed Nb library and the strategy of panning used to isolate high-affinity Nbs are shown in Figure 2.



Figure 2. Schematic representation of the construction of phage-displayed nanobody library and the strategy of panning for the high-affinity Nbs isolation.

Recently, a novel high-throughput screening strategy, "off-rate selection," was proposed by Pírez-schirme et al. [58]. Nb specific to cyanotoxin microcystin-LR (995 Da) was obtained by promoting the isolation with the slowest *koff* for the immobilized hapten and exhibited higher sensitivity than

Nbs selected by other two regular panning methods. This approach provides an alternative to improve the probability of positive anti-hapten Nb selection.

3. Performance of Nbs in the Field of Chemical Analysis

3.1. Affinity (Sensitivity)

Affinity is the binding strength between an antigenic determinant and anantigen-binding site. Early studies mostly described low-affinity Nbs [59,60]. These results are not surprising, because the overall anti-hapten sdAbs response found in Camelide are observed to be modest. It may result from only a single variable domain of sdAb contributing to form the antigen binding site, which limits the contact between the antibody and antigen, whereas conventional antibodies appear to have a relatively two-dimensional shape [61]. Table 2 provides Nb assay sensitivities (represented by 50% inhibitory concentration, IC₅₀) and, when possible, provides mAb or pAb-based assay sensitivities to the equivalent mycotoxin. The AFB_1 and OTA detection were conducted by the method of indirect competitive enzyme-linked immunosorbent assay (icELISA) based on either Nb or conventional antibodies. As for 15-acetyl-DON determination, Nb-based competitive fluorescence polarization assay was developed and compared the sensitivity with mAb/pAb-based direct competitive ELISA (dcELISA). It was noted that Nbs have similar sensitivity in comparison with the conventional antibodies. In some cases, such as Nb against 2,2',4,4'-tetrabromodiphenyl ether (BDE-47, molecular weight 485.79 Da) [62], the sensitivity of Nb-based icELISA (IC₅₀ = 1.4 ng/mL) is even higher than that of mAb ($IC_{50} = 1.75 \text{ ng/mL}$) [63]. In another study, Nb against tetrabromobisphenol-A (TBBPA, molecular weight 543.9 Da) [57], also exhibits higher sensitivity (IC₅₀ = 0.4 ng/mL) by means of the icELISA method, compared with that of mAb ($IC_{50} = 3.87 \text{ ng/mL}$) [64] and pAb $(IC_{50} = 0.87 \text{ ng/mL})$ [65].

Mycotoxin	Antibody	Assay Format	Sensitivity, IC ₅₀
AFB ₁	Nb	icELISA [55]	0.754 ng/mL
	mAb	icELISA [17]	0.0002 ng/mL
	pAb	icELISA [66]	2.0 ng/mL
OTA	Nb	icELISA [67]	0.64 ng/mL
	mAb	icELISA [68]	0.058 ng/mL
	pAb	icELISA [69]	5.0 ng/mL
15-acetyl-DON	Nb	fluorescence polarization [59]	419 ng/mL
	mAb	dcELISA [70]	1000 ng/mL
	pAb	dcELISA [71]	1.9 ng/mL

Table 2. Nbs to mycotoxins to date and comparison of sensitivities to monoclonal antibodies and polyclonal antibodies.

3.2. Selectivity

Nbs show high selectivity to macromolecular targets, especially as enzyme inhibitors due to their unique protruding CDR3, and show better selectivity toward closely related enzymes than conventional antibodies whose antigen-binding site is typically concave or flat and more likely to recognize the surface of the targets [72]. As for small haptens, it seems that no obvious difference between Nbs and conventional antibodies in terms of specificity, especially in recognizing antigens with structural analogues. The cross-reactivity of immunoassays is of usual occurrence, which is a common issue of both conventional antibodies and recombinant antibodies [60]. It was reported that the chance of selecting specific antibodies secreted from B cells mostly depends on the particular design of the immunized antigen, immunized dose, the ratio of the hapten and carrier protein and the spacer-arm lengths between the hapten and carrier protein [73–76]. However, it is feasible to produce specific Nbs towards a group of chemical analytes, such as Nb specific for OTA and its analogue, OTB,

but no significant cross-reactivity for L-phenylalanine and coumarin, which have structural similarities with parts of OTA [77].

3.3. Stability

The stability of reagents is usually a key factor when developing immunoassay products. The reagents with a long shelf-life and high stability are always preferred. It has been reported that Nbs show greater resistance to heat than conventional antibodies and are thus more suitable to be applied as immunoreagents [44]. The main distinction between the Nb fragments and conventional antibody fragments is the reversibility of the thermal unfolding process, explaining its retained functionality after exposure to high temperatures [78]. Nbs have also been observed with high stability to organic solvents, which is beneficial for developing immunoassays to be used for the analysis of lipophilic analytes, such as mycotoxins and pesticides. Since the water-miscible solvents are commonly used for extraction of analytes prior to analysis, less dilution of the extracts is needed if the matrix effects for solvent extracts are negligible, enhancing the assay's sensitivity. He et al. [55] estimated the solvent stability of Nbs against AFB₁ with methanol, dimethylsulfoxide (DMSO), dimethylformamide (DMF), acetone, and acetonitrile. By comparison of the performance of mAb against AFB₁, the results showed that Nbs displayed a higher stability to methanol, acetone, and acetonitrile than that of mAb.

3.4. Solubility

The single domain nature of Nbs leads them to facile single gene clone and easy to be expressed in soluble form in various expression systems, including bacteria, fungi, mammalian cell lines, plants, and insects [79]. The Nbs are naturally soluble because of the substitutions of hydrophobic by hydrophilic residues in the FR2. In addition, the absence of a linker peptide prone to aggregation or susceptible to proteolysis (e.g., as observed for scFv), promotes the solubility of Nb, even at concentrations above 10 mg/mL [80]. The majority of Nbs has been produced in *E. coli*, either periplasmically or cytoplasmically. Periplasmic Nbs production yields typically range from 1 to 10 mg/L while cytoplasmic Nbs production of soluble Nbs in shake flask cultures without immunizing animals, which is needed in mAbs and pAbs reproduction from host animals. The high solubility and high accumulation levels of Nbs make them ideal immunoreagents for the large-scale production in industry.

An overview of the Nbs presented within this review together with their targets, functional types, involved immunoassays and their performance are tabulated in Table 3.

ce of Nbs used for mycotoxin detection.					
Specificity (CR)	Thermal Stability	Solvent/Matrix Stability			
not tested	95 °C for 5min, retained>45% of binding activity; 90 °C for 75 min, retained>25% of binding activity	dilution factor ×2.5 eliminate the matrix interference			
1% with OTB, ZEN, DON, and AFB1	not tested	Susceptive to methanol, ionic strength and low pH (<6.0)			

Table 3. Some applications and the performance of

Targets (MW).	NB TYPE	Assay Formats, (Ref.)	Sensitivity, IC ₅₀ ; Linear Range, IC ₂₀ –IC ₈₀ ; LOD, IC ₁₀	Specificity (CR)	Thermal Stability	Solvent/Matrix Stability
OTA (403.813)	Anti-hapten Nb	icELISA [67]	IC ₅₀ : 0.64 ng/mL; LR: 0.27–1.47 ng/mL; LOD: 0.16 ng/mL	not tested	95 °C for 5min, retained>45% of binding activity; 90 °C for 75 min, retained>25% of binding activity	dilution factor ×2.5 eliminate the matrix interference
	Nb-AP	dcELISA [81]	IC ₅₀ : 0.13 ng/mL; LR: 0.06–0.43 ng/mL; LOD: 0.04 ng/mL,	CR: <0.1% with OTB, ZEN, DON, and AFB_1	not tested	Susceptive to methanol, ionic strength and low pH (\leq 6.0)
	phage-Nb	Nb-IPCR [82]	LR: 10^{-5} –1.0 ng/mL; LOD: 3.7 × 10^{-6} ng/mL;	CR: 3.5% with OTB; <0.1% with FB ₁ , ZEN, DON, and AFB ₁	not tested	not tested
	phage-AI-Nb	AI-Nb-IPCR [47]	LR: 0.01–10 ng/mL; LOD: 0.004 ng/mL	CR: <0.1% with FB ₁ , ZEN, DON, and AFB ₁	not tested	susceptive to methanol over 10%
AFB ₁ (312.27)	Anti-hapten Nb	icELISA [55]	IC ₅₀ : 0.754 ng/mL; LR: 0.117–5.676 ng/mL; LOD: 0.05 ng/mL	CR: <10% with AFB ₂ , AFG ₁ , and AFG ₂	85 °C for 60 min, retained>40% of binding activity	stable in 40% methanol and 40% acetone
	phage-AI-Nb	AI-Nb-IPCR [83]	LOD: 0.02 ng/mL	CR: >50% with AFB ₂ and AFG ₁ ; 13.5% with AFG ₂	not tested	susceptive to methanol over 10%
	AI-Nb	AI-ELISA [34]	IC ₅₀ : 0.16 ng/mL; LR: 0.09–0.82 ng/mL; LOD: 0.06 ng/mL	CR: >50% with AFB ₂ and AFG ₁ ; >30% with AFG ₂ and AFM ₁	80 °C for 60 min, retained >20% of bindingactivity	stable in methanol below 40%
	AI-Nb	TRFICA [84]	IC ₅₀ : 0.46 ng/mL; LR: 0.13–4.54 ng/mL; LOD: 0.05 ng/mL	CR: >50% with AFB ₂ and AFG ₁ ; 31.1% with AFG ₂ ; 19.4% with AFM ₁	not tested	not tested
ZEN (318.37)	phage-AI-Nb	AI-Nb-IPCR [46]	LR: 0.01–10 ng/mL; LOD: 0.0065 ng/mL	CR: <0.1% with AFB ₁ , DON, and OTA	not tested	susceptive to methanol over 5%
	AI-Nb	TRFICA [84]	IC ₅₀ : 0.86 ng/mL; LR: 0.20–2.77 ng/mL; LOD: 0.07 ng/mL;	CR: 78.1% with $\beta\text{-zearalenol}$	not tested	not tested
15-acetyl-DON (338.35)	Anti-hapten Nb	fluorescence polarization [59]	IC ₅₀ : 419.5 ng/mL (monomer) 167.175 ng/mL (pentamer)	CR: <0.1% with neosolaniol, diacetoxyscirpenol, and T-2 toxin	not tested	not tested
DON (296.32)	AI-Nb	AI-ELISA [48]	IC ₅₀ : 8.77 ng/mL; LR: 2.18–62.25 ng/mL; LOD: 1.16 ng/mL	CR: <0.1% with FB ₁ , ZEN, AFB ₁ , and OTA	95 °C for 5 min, retained>60% of binding activity	not tested
	phage-AI-Nb	AI-PELISA [85]	IC ₅₀ : 24.49 ng/mL; LR: 9.51–180.15 ng/mL	CR: <0.1%with FB ₁ , ZEN, AFB ₁ , and OTA	95 °C for 5 min, retained>55% of binding activity	stable in ionic strength below 50 mM and pH between 5.0 and 8.0

Targets (MW).	NB TYPE	Assay Formats, (Ref.)	Sensitivity, IC ₅₀ ; Linear Range, IC ₂₀ –IC ₈₀ ; LOD, IC ₁₀	Specificity (CR)	Thermal Stability	Solvent/Matrix Stability
FB ₁ (721.84)	AI-Nb	AI-ELISA [49]	IC ₅₀ : 0.95 ng/mL; LR: 0.27–5.92 ng/mL; LOD: 0.15 ng/mL	CR: 4.93% with FB ₂	not tested	stable in ionic strength below 50 mM and pH between 5.0 and 8.0
CIT (250.25)	AI-Nb	AI-ELISA [50]	IC ₅₀ : 44.6 ng/mL; LR: 5.0–300.0 ng/mL;	CR: <0.1% with AFB ₁ , ZEN, DON, and OTA	not tested	susceptive to methanol over 20%
	phage-AI-Nb	AI-PELISA [45]	IC ₅₀ : 10.9 ng/mL; LR: 2.5–100.0 ng/mL	CR: <0.1% with AFB ₁ , ZEN, DON, and OTA	not tested	stable in methanol below 25%, ionic strength below 55 mM and pH between 5.4 and 9.0

Abbreviations used: (1) MW, molecular weight; (2) IC₅₀, IC₂₀–IC₈₀, IC₁₀, concentrations resulting in 50, 20-80 or 10% decrease in maximum signal; (3) LR, linear range; (4) LOD, limit of detection; (5) CR, cross-reactivity; (6) ELISA, enzyme linked immunosorbent assay; (7) ic-ELISA, indirect competitive ELISA; (8) dc-ELISA, direct competitive ELISA; (9) Nb-IPCR, Nb-based phage displayed immune PCR; (10) AI-Nb-IPCR, anti-idiotypic Nb based phage displayed immune PCR; (11) AI-ELISA, anti-idiotypic Nb based ELISA; (12) AI-PELISA, anti-idiotypic Nb-based phage ELISA; (13) TRFICA, time-resolved fluorescence immunochromatographic assay; (14) Nb-AP, Nb-alkaline phosphatase fusion protein; (15) phage-Nb, phage-displayed Nb; (16) phage-AI-Nb, phage-displayed anti-idiotypic Nb. * units converted to ng/mL from initial publication.

4. Application of Nbs in the Area of Mycotoxin Detection

In spite of the unfavorable camelid immune response against small haptens, high-affinity Nbs to a diverse of targets have been generated and utilized to develop highly sensitive immunoassays for mycotoxin detection. Those methods, including Nb-based icELISA, one-step dcELISA based on Nb-alkaline phosphatase (AP) fusion protein, AI-Nb-based ELISA (AI-ELISA) and phage-Nb combined with immune-PCR. Some other exquisite designed immunoassays, such as Nb-based fluorescence polarization immunoassay and immunochromatographic assay, were also proved to be feasible, which widens the Nb application in mycotoxin detection in the area of food quality control.

4.1. Competitive Enzyme-Linked Immunosorbent Aassay

There are two normal types of competitive ELISA, indirect and direct competitive assay, applied in small chemical determination. Figure 3 displays the main forms of Nb-based competitive ELISA.



(2) Competitive direct ELISA format, C and D

Figure 3. Schematic diagram of Nb-based competitive ELISA formats. The direct assays do not use second antibody and have one less step than the indirect assays. **(A)** Nb-based icELISA; **(B)** anti-idiotypic Nb-based icELISA, AI-ELISA; **(C)** Nb-based dcELISA; **(D)** anti-idiotypic Nb-based dcELISA, AI-dcELISA.

Anti-hapten Nbs can be used as affinity reagents in immune recognition and to construct a Nb-based icELISA (Figure 3A). Liu et al. [67] raised positive Nb against OTA and developed a sensitive Nb-based icELISA. The assay showed high resistance to cereal matrix interference and was applied to OTA-contaminated cereal samples determination. He et al. [55] developed an icELISA based on Nb for the analysis of AFB₁, in which Nb also showed high stability under organic solvents and high temperature. The recovery of the assay was acceptable for the quantitative analysis of AFB₁ in agro-product samples without dilution of extracted samples.

With the feasibility of genetic manipulation, Nb can be easily combined with alkaline phosphatase (Nb-AP) as fusion protein for direct colorimetric detection and construct the dcELISA (Figure 3C), which can avoid the chemical conjugation of a second enzyme and reduce the reaction step without altering sensitivity. Liu et al. [81] developed dcELISA based on a Nb-AP fusion protein for monitoring OTA contamination in cereal. The IC₅₀ of the assay was 0.13 ng/mL, which increased 5-fold compared to Nb-based icELISA. Moreover, the detection limit and linear range were 0.04 ng/mL and 0.06-0.43 ng/mL, respectively. Shu et al. [49] developed a green and rapid one-step competitive enzyme

immunoassay based on an AI-Nb and AP as fusion protein for FB₁ detection in cereal (Figure 3D). The sensitivity of one-step dcELISA was 9-fold higher than two-step icELISA when the same isolated AI-Nb was involved in two types of ELISA.

4.2. Anti-Idiotypic Nb-based ELISA

Anti-idiotypic antibodies are novel antigens and generated according to primary antibodies, giving rise to humoral immune responses in syngeneic or xenogeneic systems [86]. The second antibody targets the antigenic determinants of the primary antibody, which can replace the conventional hapten-protein conjugates of small molecules, serving the same function in the competitive immunoassay. The Nbs against primary antibodies can be used as anti-idiotypic antibodies and work as non-toxic candidates for competing antigens in immunoassays [87] (Figure 3B).

Wang et al. [34] isolated AI-Nb from anti-aflatoxin mAb immunized alpaca library and then applied to immunoassay (AI-ELISA) toward aflatoxin detection as a coating antigen. The best immunoassay developed with the AI-Nb showed an IC_{50} value of 0.16 ng/mL toward AFB₁ and cross-reactivity toward aflatoxin B₂, G₁, G₂, and M₁ of 90.4%, 54.4%, 37.7%, and 37.4%, respectively. The AI-ELISA was successfully applied to the aflatoxin analysis of peanuts, corn, and rice. The final assay has a good sensitivity and can be used to detect AFB₁ in agro-products. The newly developed assay avoids the use of an aflatoxin standard in the synthesis of a hapten-protein conjugate as a coating antigen, which was replaced by an AI-Nb.

Xu et al. [45,50] developed a simple and straightforward strategy to isolate AI-Nb as a CIT mimotope from a naive phage display Nb library. The obtained AI-Nb was successfully used as a competing antigen in establishing AI-ELISA and phage ELISA (AI-PELISA) for detection of CIT. The AI-ELISA exhibited dynamic linear detection of CIT from 5.0 to 300.0 ng/mL, with an IC₅₀ value of 44.6 ng/mL, which was twice as good as that of icELISA (IC₅₀ = 96.2 ng/mL) with CIT-OVA as a coating antigen. The AI-PELISA performed dynamic linear range from 2.5 ng/mL to 100.0 ng/mL with an IC₅₀ value of 10.9 ng/mL, which was 9 times more sensitive than conventional CIT-OVA conjugates-based icELISA.

Qiu et al. [48] selected AI-Nb, N-28, which bind to anti-DON mAb from a naive phage display library and developed an AI-ELISA for sensitive detection of DON in cereals and feedstuffs. The IC₅₀ value of the immunoassay developed was 8.77 ng/mL, which was 18-fold more sensitive than the conventional coating antigen (DON-BSA) based immunoassay. In addition, the interaction mechanism between isolated AI-Nb and anti-DON mAb was characterized and validated by molecular modeling and alanine-scanning mutagenesis. The results showed that hydrogen bond and hydrophobic interaction may contribute to their binding. Then, Qiu et al. [85] conducted site-saturation mutagenesis at five residues of N-28 and constructed a mutant phage display library. The AI-PELISA based on the best mutant N-28-T102Y was established and the IC₅₀ value of the assay (IC₅₀ = 24.49 ng/mL) was 3.2 times more sensitive than the wide-type N-28-based AI-PELISA (IC₅₀ = 78.82 ng/mL). The linear range of the assay was 9.51–180.15 ng/mL. The established AI-PELISA was applied into DON detection in wheat, corn and feedstuff samples with good recovery values and the results of the real sample detection by AI-PELISA were in agreement with commercial ELISA kit.

In another work, Shu et al. [49] made use of an anti-FB₁ mAb as the target for panning FB₁ substitute from a naive alpaca Nb phage display library. The isolated AI-Nb was subjected to an AI-ELISA for the detection of FB₁ contaminated in cereals and feedstuffs. The developed assay showed an IC₅₀ value of 0.95 ng/mL, with a limit of detection of 0.15 ng/mL, linear range of 0.27–5.92 ng/mL, and low cross-reactivity toward FB₂ (4.93%). The sensitivity of AI-ELISA was enhanced approximately 20-fold compared with that of the chemosynthetic FB₁-BSA conjugates-based ELISA (IC₅₀ = 21.14 ng/mL). The established AI-ELISA was validated to be suitable for monitoring the total fumonisin concentration under the current regulatory limits of fumonisins in most countries.

4.3. Phage Display Mediated Immuno-Polymerase Chain Reaction

Phage display mediated immune-PCR (PD-IPCR) has been proved to be an ultrasensitive combined method using immunoassay and PCR for determination of toxins and markers. The phage-Nb and phage DNA can directly act as the detection antibody and PCR template, respectively.

There are two forms of phage-borne Nb used to construct real-time IPCR (RT-IPCR) assay, one of which is antigen-specific Nb phage-based RT-IPCR (Nb-IPCR) and the other is AI-Nb phage-based RT-IPCR (AI-Nb-IPCR). The former phage-Nb is able to recognize immobilized antigen and compete with free analytes while the latter AI-Nb phage (phage-AI-Nb) against anti-hapten mAb can be used as surrogate for competitive antigen. Both of these two types of phage-displayed Nbs (phage-Nb and phage-AI-Nb) can be used as signal-amplification carrier with specific recombinant DNA to develop an environmentally friendly IPCR for chemical contaminant detection in agriculture and environment. Figure 4 shows the schematic diagram of the PD-IPCR assay format.



Figure 4. Schematic diagram of the real-time phage display mediated immune-PCR assay. (Reprinted with permission from reference [83], Copyright (2014), American Chemical Society.)

Liu et al. [82] generated OTA-specific Nb from an immunized alpaca Nb library and then applied phage-Nb to IPCR for ultrasensitive detection of OTA in cereal samples. The detection limit of the Nb-IPCR was 3.7 pg/L (LOD = 3.7×10^{-6} ng/mL), and with a wide linear range of 0.01–1000 pg/mL (LR: 10^{-5} –1.0 ng/mL). The validation results indicated the reliability of Nb-IPCR in the detection of OTA in cereal samples. This method is the most sensitive assay reported to date for the detection of OTA and indicates that Nb-IPCR has great potential in the ultrasensitive detection of mycotoxins and other toxic small molecular compounds.

In another study, Ji et al. [47] developed a sensitive and non-toxin quantitative AI-Nb-IPCR for OTA based on phage-AI-Nb against anti-OTA mAb. The limit of detection of assay was 4.17 pg/mL, which exhibited a 9-fold improvement over AI-PELISA (LOD = 300 pg/mL). Moreover, a wide linear range (LR: 10–10,000 pg/mL) was available. The developed method was successfully validated with OTA contaminated agro-products.

Lei et al. [83] developed AI-Nb-IPCR assay for the accurately quantitative detection of major aflatoxins (AFB₁, B₂, G₁ and G₂) in agro-products based on phage-AI-Nb against anti-aflatoxins mAb. The limit of detection of the established assay was 0.02 ng/mL, which was 4-fold improvement over AI-PELISA (LOD = 0.08 ng/mL) for AFB₁ detection. The developed method was successfully validated with the samples of corn, rice, peanut, and feedstuff. Then, Wang et al. [44] also developed AI-Nb-IPCR for ultrasensitive determination of ZEN in cereals. The limit of detection of the assay was 6.5 pg/mL,

which was 12 times improved than that of exclusive AI-PELISA (LOD = 0.25 ng/mL), with a linear range of 0.01-100 ng/mL. The proposed method was successfully applied for the determination of ZEN in cereal samples.

4.4. Nb-Based Fluorescence Polarization Immunoassay

Doyle et al. [59] isolated Nb against 15-acetyl-DON from a hyper-immunized llama library, followed by soluble expression of monomeric and pentameric Nb for application in fluorescence polarization immunoassay. DON was conjugated with fluorescein and used as tracer. Serial dilutions of Nb were mixed with tracer to react in the dark and followed by fluorescence polarization analysis with excitation and emission wavelengths of 488 and 530 nm, respectively. A standard curve of competition between DON-fluorescein and free DON was prepared using a competitive fluorescence polarization assay. The assay inhibition with 15-acetyl-DON determined an IC₅₀ value of 1.24 μ M (419.5 ng/mL) and 0.50 μ M (167.175 ng/mL) for monomer and pentamer Nb, respectively.

4.5. Nb-based Lateral-Flow Immunoassays

Antibody-based immunochromatographic assay is widely used in rapid detection of mycotoxins in agro-products. Tang et al. [84] developed time-resolved fluorescence immunochromatographic assay (TRFICA) using two AI-Nbs for rapid, quantitative, and simultaneous detection of AFB₁ and ZEN in maize and its products. A novel lanthanides Eu/Tb(III) nanosphere with enhanced fluorescence was conjugated to two anti-hapten mAbs—one is against AFB₁ and the other is against ZEN—and used as detectors. Two AI-Nbs, AFB₁ surrogate and ZEN surrogate, serving as coating antigen and immobilized on two separating test lines as capture antigens. The pattern of competitive time-resolved strip methods (AI-Nb-TRFICA) was established and the standard carve was built. The IC₅₀ value of AI-Nb-TRFICA was 0.46 and 0.86 ng/mL for AFB₁ and ZEN, respectively. The simultaneous detection of dual mycotoxins by AI-Nb-TRFICA was established, which provided a quantitative relationship ranging from 0.13 to 4.54 ng/mL for AFB₁ and 0.20 to 2.77 ng/mL for ZEN, with a detection limit of 0.05 and 0.07 ng/mL in the buffer solution, respectively. The assay showed good recoveries in samples and was applied to detect dual mycotoxins in maize samples with satisfied results. Figure 5 displays the schematic diagram of the AI-Nb-based TRFICA assay format.



Figure 5. Schematic illustration of AI-Nb-based time-resolved fluorescence immunochromatographic assay for dual mycotoxins. (Reprinted with permission from reference [84], Copyright (2017), American Chemical Society.)

This research is the first report about time-resolved strip method based on AI-Nb for dual mycotoxins detection, which provides evidence that the Nb reagent possesses a broad potential for a wide range of application in the field of food quality control in the future.

4.6. Nb-Based Immunoaffinity Method

Sample pretreatment and purification are usually needed previous to pollutants analysis. Antibody-based immunoaffinity method is a common tool used for samples clean-up. Nb has been proved to possess high degeneration resistance that would enable the purification reagent reuse and further replace conventional antibodies. Xiong et al. [88] applied anti-AFB₁ Nb into magneticbead-based immunoaffinity extraction (M-IAE) method for purification of AFB₁ from corn samples. Magnetic beads carrying poly (acrylic acid) brushes (MB@PAA) were fabricated as an "Nb container" for improving AFB₁ adsorption capacity. The developed MB@PAA@Nb could be reused at least 10 times, without obvious loss of the capture efficiency for AFB₁. This proposed MB@PAA@Nbs-based immunoaffinity extraction method is a highly promising, novel sample pre-treatment platform for AFB₁ as well as other mycotoxins.

5. Challenge of Nb Application in Analytical Chemistry

Nbs have many favorable properties compared with conventional and other recombinant antibodies. The smaller number of CDRs does not seem to threaten the ability of Nbs to bind macromolecular antigens with high affinity, such as AI-Nb binding to mAb is easy to be obtained. However, the lack of a light chain possibly hampers their ability to bind small haptens that conventional antibodies typically recognize by the interplay of the VH and VL domains and six CDRs [89]. This is the case for the isolation of Nbs against small molecules, like mycotoxins.

Recently, only three anti-mycotoxin Nbs are generated. Failures in the generation of Nbs that recognize small haptens in solution are common [43]. It seems difficult to get out of dilemma even sophisticated screening routine operated [60]. The reason may be the camelid immune response against small molecules is predominated by conventional IgG1, with lower ratios of single domain subclasses of IgG2 and IgG3 [54]. Additional structural analysis of small molecule-Nb complexes would be helpful in understanding the interaction and leading a more rational design of immunizing and selecting haptens. However, there is very little information on the way in which Nbs bind small molecules. So far, only three crystal structures of hapten-Nb complexes have been revealed, corresponding to the azo dyes Reactive Red 1 (733 Da) [90], Reactive Red 6 (717 Da) [91], and the chemotherapy agent methotrexate (454 Da) [92]. Among these three complexes, all Nbs can provide cavities to accommodate haptens like the VH-VL crevice, even in the absence of a VL. The three CDRs usually contribute to form the groove except methotrexate, which is bound in a non-canonical way, being deeply buried in a "tunnel" roofed by CDR1 and a loop of FR3 that the authors identified as "CDR4". As a result, more efforts are needed to generate Nbs against small haptens and to further determine the structure characteristic and molecular recognition mechanism. Owing to their small size and ease of expression, Nbs are usually easy to crystallize or operate nuclear magnetic resonance (NMR) determination for structure study. It suggests that more Nb structures may be elucidated soon. Those emerging structures will help to identify the possible mechanisms and patterns of the interaction between Nbs and their targets, especially as to Nbs recognizing small haptens. Understanding the structure of the combining pocket will help to adjust Nbs characteristics through in vitro antibody maturation techniques, and in addition, to guide artificial antigens design rationally.

Additionally, Nb-based immunodetection formats are still limited and need to be further explored for mycotoxin detection, such as biosensors and portable devices [93,94]. It has been reported that the anti-hapten Nbs can be employed as recognition element and incorporated into two biosensor formats. The Nb can be coated on the surface of an electrode for electrochemical impedance detection [62], and coated onto a PDMS membrane for a lab-on-a-chip sensor [95]. These studies demonstrate the potential of developing novel and reliable assays for mycotoxin detection by Nb technology.

6. Conclusions

Ever since their discovery, Nbs have emerged as powerful antigen binders and have become a very promising tool for the monitoring of food quality. In this review, we summarized the recent advancements of Nb-based immunoassays for mycotoxin detection, including indirect competitive ELISA, direct competitive ELISA, phage-displayed Nb-mediated immuno-PCR, fluorescence polarization immunoassay, immunochromatographic assay, and immunoaffinity methods. Lots of studies have demonstrated Nb technology to be a powerful and potential technique for advanced applications in mycotoxin detection. As already mentioned, Nbs have striking properties such as high solubility, ease of genetic manipulation, and high thermal and conformational stability. However, it is still very challenging to select desirable Nbs binding to the small molecule epitopes from the phage display Nb repertoire, hampering the more extensive applications of Nbs in mycotoxin detection. Efforts should be made as follows to cope with the difficulties: (1) improve the success rate of Nb generation. It is necessary to ensure the immunogenicity of the designed hapten-protein conjugate, the capacity and diversity of the constructed Nb repertoire, as well as the optimization of the screening strategy; (2) the elucidation of the structure of Nbs and their complex with small molecules is a promising way to modify the properties of Nbs and to provide clues to the rational design of immunogen for the anti-hapten Nbs generation; and (3) innovative immunodetection formats should be further exploited to improve the performance of Nb-based assay and to expand the applicability of Nbs for the mycotoxin detection.

Author Contributions: Y.Y. and P.L. conceived and designed the work; Q.Z. co-designed the work and revised the manuscript; T.W. and Y.N. collected the references; J.Z. and R.H. analyzed the data; T.H. wrote the paper.

Acknowledgments: This work was supported by funds from the National Natural Sciences Foundation of China (Grant Number: 21575155, 31471650), and the Hundred Talent Program by Chinese Academy of Sciences.

Conflicts of Interest: All authors declare no conflict of interest.

References

- Stoev, S.D. Food safety and increasing hazard of mycotoxin occurrence in foods and feeds. *Crit. Rev. Food Sci.* 2013, 53, 887–901. [CrossRef] [PubMed]
- Hussein, H.S.; Brasel, J.M. Toxicity, metabolism, and impact of mycotoxins on humans and animals. *Toxicology* 2001, 167, 101–134. [CrossRef]
- Sweeney, M.J.; Dobson, A.D.W. Mycotoxin production by Aspergillus, Fusarium and Penicillium species. *Int. J. Food Microbiol.* 1998, 43, 141–158. [CrossRef]
- Shuib, N.S.; Makahleh, A.; Salhimi, S.M.; Saad, B. Determination of aflatoxin M1 in milk and dairy products using high performance liquid chromatography-fluorescence with post column photochemical derivatization. *J. Chromatogr. A* 2017, 1510, 51–56. [CrossRef] [PubMed]
- Kim, D.H.; Hong, S.Y.; Kang, J.W.; Cho, S.M.; Lee, K.R.; An, T.K.; Lee, C.; Chung, S.H. Simultaneous determination of multi-mycotoxins in cereal grains collected from South Korea by LC/MS/MS. *Toxins* 2017, 9, 147. [CrossRef] [PubMed]
- 6. Urusov, A.E.; Zherdev, A.V.; Dzantiev, B.B. Immunochemical methods of mycotoxin analysis (review). *Appl. Biochem. Microbiol.* **2010**, *46*, 253–266. [CrossRef]
- Tripathi, P.; Upadhyay, N.; Nara, S. Recent advancements in lateral flow immunoassays: A journey for toxin detection in food. *Crit. Rev. Food Sci.* 2017. [CrossRef] [PubMed]
- 8. Tang, X.Q.; Li, X.; Li, P.W.; Zhang, Q.; Li, R.; Zhang, W.; Ding, X.X.; Lei, J.W.; Zhang, Z.W. Development and application of an immunoaffinity column enzyme immunoassay for mycotoxin zearalenone in complicated samples. *PLoS ONE* **2014**, *9*, e85606. [CrossRef] [PubMed]
- Xue, S.; Li, H.P.; Zhang, J.B.; Liu, J.L.; Hu, Z.Q.; Gong, A.D.; Huang, T.; Liao, Y.C. Chicken single-chain antibody fused to alkaline phosphatase detects aspergillus pathogens and their presence in natural samples by direct sandwich enzyme-linked immunosorbent assay. *Anal. Chem.* 2013, *85*, 10992–10999. [CrossRef] [PubMed]

- Wang, T.; Li, P.W.; Zhang, Q.; Zhang, W.; Zhang, Z.W.; Wang, T.; He, T. Determination of aspergillus pathogens in agricultural products by a specific nanobody-polyclonal antibody sandwich ELISA. *Sci. Rep.* 2017, *7*, 4348. [CrossRef] [PubMed]
- Li, P.W.; Zhang, Q.; Zhang, W. Immunoassays for aflatoxins. *Trac-Trends Anal. Chem.* 2009, 28, 1115–1126. [CrossRef]
- 12. EU No. 165/2010. Commission Regulation (EU) No. 165/2010 of 26 February 2010, Amending Regulation (EC) No. 1881/2006, Setting Maximum Levels for Certain Contaminants in Foodstuffs as Regards Aflatoxins (Text with EEA Relevance); European Commission: Brussels, Belgium, 2010.
- EC No. 1126/2007. Commission Regulation (EC) No. 1126/2007 of 28 September 2007, Amending Regulation (EC) No. 1881/2006, Setting Maximum Levels for Certain Contaminants in Foodstuffs as Regards Fusarium Toxins in Maize and Maize Products (Text with EEA Relevance); European Commission: Brussels, Belgium, 2007.
- EU No 212/2014. Commission Regulation (EU) No. 212/2014of 6 March 2014, Amending Regulation (EC) No. 1881/2006 as Regards Maximum Levels of the Contaminant Citrinin in Food Supplements Based on Rice Fermented with Red Yeast Monascuspurpureus (Text with EEA Relevance); European Commission: Brussels, Belgium, 2014.
- Shiu, C.M.; Wang, J.J.; Yu, F.Y. Sensitive enzyme-linked immunosorbent assay and rapid one-step immunochromatographic strip for fumonisin B1 in grain-based food and feed samples. *J. Sci. Food Agric.* 2010, *90*, 1020–1026. [CrossRef] [PubMed]
- Zhuang, Z.H.; Que, S.J.; Gao, Y.M.; Yuan, J.; Ye, Z.; Du, M.; Lin, G.M.; Liu, L.C.; Wang, S.H. Artificial antigen synthesis and the development of polyclonal antibody-based immunoassay for citreoviridin determination. *World J. Microbiol. Biotechnol.* 2014, 30, 343–349. [CrossRef] [PubMed]
- 17. Zhang, D.H.; Li, P.W.; Zhang, Q.; Zhang, W.; Huang, Y.L.; Ding, X.X.; Jiang, J. Production of ultrasensitive generic monoclonal antibodies against major aflatoxins using a modified two-step screening procedure. *Anal. Chim. Acta* **2009**, *636*, 63–69. [CrossRef] [PubMed]
- KalayuYirga, S.; Ling, S.; Yang, Y.; Yuan, J.; Wang, S. The preparation and identification of a monoclonal antibody against citrinin and the development of detection via indirect competitive ELISA. *Toxins* 2017, 9, 110. [CrossRef] [PubMed]
- 19. Li, X.; Li, P.W.; Lei, J.W.; Zhang, W.; Li, C.M. A simple strategy to obtain ultra-sensitive single-chain fragment variable antibodies for aflatoxin detection. *Rsc Adv.* **2013**, *3*, 22367–22372. [CrossRef]
- 20. Arola, H.O.; Tullila, A.; Nathanail, A.V.; Nevanen, T.K. A simple and specific noncompetitive ELISA method for HT-2 toxin detection. *Toxins* **2017**, *9*, 145. [CrossRef] [PubMed]
- 21. Hamers-Casterman, C.; Atarhouch, T.; Muyldermans, S.; Robinson, G.; Hammers, C.; Songa, E.B.; Bendahman, N.; Hammers, R. Naturally occurring antibodies devoid of light chains. *Nature* **1993**, *363*, 446–448. [CrossRef] [PubMed]
- 22. Narciso, J.E.T.; Uy, I.D.C.; Cabang, A.B.; Chavez, J.F.C.; Pablo, J.L.B.; Padilla-Concepcion, G.P.; Padlan, E.A. Analysis of the antibody structure based on high-resolution crystallographic studies. *New Biotechnol.* **2011**, *28*, 435–447. [CrossRef] [PubMed]
- 23. Arbabi Ghahroudi, M.; Desmyter, A.; Wyns, L.; Hamers, R.; Muyldermans, S. Selection and identification of single domain antibody fragments from camel heavy-chain antibodies. *FEBS Lett.* **1997**, 414, 521–526. [CrossRef]
- Cortez-Retamozo, V.; Lauwereys, M.; Gh, G.H.; Gobert, M.; Conrath, K.; Muyldermans, S.; Baetselier, P.D.; Revets, H. Efficient tumor targeting by single-domain antibody fragments of camels. *Int. J. Cancer* 2002, *98*, 456–462. [CrossRef] [PubMed]
- 25. Hudson, P.J.; Souriau, C. Engineered antibodies. Nat. Med. 2003, 9, 129–134. [CrossRef] [PubMed]
- 26. Yau, K.Y.; Lee, H.; Hall, J.C. Emerging trends in the synthesis and improvement of hapten-specific recombinant antibodies. *Biotechnol. Adv.* 2003, 21, 599–637. [CrossRef]
- 27. Le Gall, F.; Reusch, U.; Moldenhauer, G.; Little, M.; Kipriyanov, S.M. Immunosuppressive properties of anti-CD3 single-chain Fv and diabody. *J. Immunol. Methods* **2004**, *285*, 111–127. [CrossRef] [PubMed]
- 28. Vu, K.B.; Ghahroudi, M.A.; Wyns, L.; Muyldermans, S. Comparison of llama VH sequences from conventional and heavy chain antibodies. *Mol. Immunol.* **1997**, *34*, 1121–1131. [CrossRef]

- Chan, P.H.; Pardon, E.; Menzer, L.; De Genst, E.; Kumita, J.R.; Christodoulou, J.; Saerens, D.; Brans, A.; Bouillenne, F.; Archer, D.B.; et al. Engineering a camelid antibody fragment that binds to the active site of human lysozyme and inhibits its conversion into amyloid fibrils. *Biochemistry* 2008, 47, 11041–11054. [CrossRef] [PubMed]
- LutjeHulsik, D.; Liu, Y.Y.; Strokappe, N.M.; Battella, S.; El Khattabi, M.; Mccoy, L.E.; Sabin, C.; Hinz, A.; Hock, M.; Macheboeuf, P.; et al. A gp41 MPER-specific llama VHH requires a hydrophobic CDR3 for neutralization but not for antigen recognition. *PLoS Pathog.* 2013, *9*, e1003202. [CrossRef]
- Muyldermans, S.; Baral, T.; Retamozzo, V.C.; De Baetselier, P.; De Genst, E.; Kinne, J.; Leonhardt, H.; Magez, S.; Nguyen, V.; Revets, H. Camelid immunoglobulins and nanobody technology. *Vet. Immunol. Immunopathol.* 2009, 128, 178–183. [CrossRef] [PubMed]
- 32. Saerens, D.; Frederix, F.; Reekmans, G.; Conrath, K.; Jans, K.; Brys, L.; Huang, L.; Bosmans, E.; Maes, G.; Borghs, G.; et al. Engineering camel single-domain antibodies and immobilization chemistry for human prostate-specific antigen sensing. *Anal. Chem.* **2005**, *77*, 7547–7555. [CrossRef] [PubMed]
- 33. Franco, E.J.; Sonneson, G.J.; Delegge, T.J.; Hofstetter, H.; Horn, J.R.; Hofstetter, O. Production and characterization of a genetically engineered anti-caffeine camelid antibody and its use in immunoaffinity chromatography. *J. Chromatogr. B* **2010**, *878*, 177–186. [CrossRef] [PubMed]
- 34. Wang, Y.; Li, P.; Majkova, Z.; Bever, C.R.S.; Kim, H.J.; Zhang, Q.; Dechant, J.E.; Gee, S.J.; Hammock, B.D.A. Chem, Isolation of alpaca anti-idiotypic heavy-chain single-domain antibody for the aflatoxin immunoassay. *Anal. Chem.* **2013**, *85*, 8298–8303. [CrossRef] [PubMed]
- 35. Conrath, K.E.; Lauwereys, M.; Wyns, L.; Muyldermans, S. Camel single-domain antibodies as modular building units in bispecific and bivalent antibody constructs. *J. Biol. Chem.* **2001**, *276*, 7346–7350. [CrossRef] [PubMed]
- De Groeve, K.; Deschacht, N.; De Koninck, C.; Caveliers, V.; Lahoutte, T.; Devoogdt, N.; Muyldermans, S.; De Baetselier, P.; Raes, G. Nanobodies as tools for in vivo imaging of specific immune cell types. *J. Nucl. Med.* 2010, 51, 782–789. [CrossRef] [PubMed]
- 37. Schoonooghe, S.; Laoui, D.; Van Ginderachter, J.A.; Devoogdt, N.; Lahoutte, T.; De Baetselier, P.; Raes, G. Novel applications of nanobodies for in vivo bio-imaging of inflamed tissues in inflammatory diseases and cancer. *Immunobiology* **2012**, *217*, 1266–1272. [CrossRef] [PubMed]
- 38. Oliveira, S.; Heukers, R.; Sornkom, J.; Kok, R.J.; van Bergen EnHenegouwen, P.M. Targeting tumors with nanobodies for cancer imaging and therapy. *J. Control. Release* **2013**, *172*, 607–617. [CrossRef] [PubMed]
- Yau, K.Y.; Groves, M.A.; Li, S.; Sheedy, C.; Lee, H.; Tanha, J.; Mackenzie, C.R.; Jermutus, L.; Hall, J.C. Selection of hapten-specific single-domain antibodies from a non-immunized llama ribosome display library. *J. Immunol. Methods* 2003, 281, 161–175. [CrossRef] [PubMed]
- 40. Sheedy, C.; Yau, K.Y.; Hirama, T.; Mackenzie, C.R.; Hall, J.C. Selection, characterization, and CDR shuffling of naive llama single-domain antibodies selected against auxin and their cross-reactivity with auxinic herbicides from four chemical families. *J. Agric. Food Chem.* **2006**, *54*, 3668–3678. [CrossRef] [PubMed]
- Ladenson, R.C.; Crimmins, D.L.; Landt, Y.; Ladenson, J.H. Isolation and characterization of a thermally stable recombinant anti-caffeine heavy-chain antibody fragment. *Anal. Chem.* 2006, 78, 4501–4508. [CrossRef] [PubMed]
- 42. Alvarez-Rueda, N.; Behar, G.; Ferré, V.; Pugniere, M.; Roquet, F.; Gastinel, L.; Jacquot, C.; Aubry, J.; Baty, D.; Barbet, J. Generation of llama single-domain antibodies against methotrexate, a prototypical hapten. *Mol. Immunol.* **2007**, *44*, 1680–1690. [CrossRef] [PubMed]
- Bever, C.S.; Dong, J.X.; Vasylieva, N.; Barnych, B.; Cui, Y.; Xu, Z.L.; Hammock, B.D.; Gee, S.J. VHH antibodies: Emerging reagents for the analysis of environmental chemicals. *Anal. Bioanal. Chem.* 2016, 408, 1–18. [CrossRef] [PubMed]
- 44. Wang, Y.; Li, P.; Zhang, Q.; Hu, X.; Zhang, W. A toxin-free enzyme-linked immunosorbent assay for the analysis of aflatoxins based on a VHH surrogate standard. *Anal. Bioanal. Chem.* **2016**, *408*, 1–8. [CrossRef] [PubMed]
- 45. Xu, Y.; Xiong, L.; Li, Y.; Xiong, Y.; Tu, Z.; Fu, J.; Tang, X. Citrinin detection using phage-displayed anti-idiotypic single-domain antibody for antigen mimicry. *Food Chem.* **2015**, *177*, 97–101. [CrossRef] [PubMed]
- Wang, X.; He, Q.; Xu, Y.; Liu, X.; Shu, M.; Tu, Z.; Li, Y.; Wang, W.; Cao, D. Anti-idiotypic VHH phage display-mediated immuno-PCR for ultrasensitive determination of mycotoxin zearalenone in cereals. *Talanta* 2016, 147, 410–415. [CrossRef] [PubMed]

- 47. Ji, Y.; He, Q.; Xu, Y.; Tu, Z.; Yang, H.; Qiu, Y.; Wang, X.; Liu, Y. Phage displayed anti-idiotypic nanobody mediated immuno-PCR for sensitive and environmentally friendly detection of mycotoxin ochratoxin A. *Anal. Methods* **2016**, *8*, 7824–7831. [CrossRef]
- 48. Qiu, Y.L.; He, Q.H.; Xu, Y.; Bhunia, A.K.; Tu, Z.; Chen, B.; Liu, Y.Y. Deoxynivalenol-mimic nanobody isolated from a naïve phage display nanobody library and its application in immunoassay. *Anal. Chim. Acta* **2015**, *887*, 201–208. [CrossRef] [PubMed]
- Shu, M.; Xu, Y.; Wang, D.; Liu, X.; Li, Y.; He, Q.; Tu, Z.; Qiu, Q.; Ji, Y.; Wang, X. Anti-idiotypic nanobody: A strategy for development of sensitive and green immunoassay for fumonisin B1. *Talanta* 2015, 143, 388–393. [CrossRef] [PubMed]
- Xu, Y.; Xiong, L.; Li, Y.; Xiong, Y.; Tu, Z.; Fu, J.; Chen, B. Anti-idiotypic nanobody as citrinin mimotope from a naive alpaca heavy chain single domain antibody library. *Anal. Bioanal. Chem.* 2015, 407, 5333–5341. [CrossRef] [PubMed]
- Pardon, E.; Laeremans, T.; Triest, S.; Rasmussen, S.G.; Wohlkönig, A.; Ruf, A.; Muyldermans, S.; Hol, W.G.; Kobilka, B.K.; Steyaert, J. A general protocol for the generation of Nanobodies for structural biology. *Nat. Protoc.* 2014, *9*, 674–693. [CrossRef] [PubMed]
- 52. Desmyter, A.; Spinelli, S.; Roussel, A.; Cambillau, C. Camelid nanobodies: Killing two birds with one stone. *Curr. Opin. Struct. Biol.* **2015**, *32*, 1–8. [CrossRef] [PubMed]
- Kim, H.J.; McCoy, M.R.; Majkova, Z.; Dechant, J.E.; Gee, S.J.; Tabares-da Rosa, S.; González-Sapienza, G.G.; Hammock, B.D. Isolation of alpaca anti-hapten heavy chain single domain antibodies for development of sensitive immunoassay. *Anal. Chem.* 2011, *84*, 1165–1171. [CrossRef] [PubMed]
- 54. Tabares-da Rosa, S.; Rossotti, M.; Carleiza, C.; Carrión, F.; Pritsch, O.; Ahn, K.C.; Last, J.A.; Hammock, B.D.; González-Sapienza, G. Competitive selection from single domain antibody libraries allows isolation of high-affinity antihapten antibodies that are not favored in the llama immune response. *Anal. Chem.* 2011, *83*, 7213–7220. [CrossRef] [PubMed]
- He, T.; Wang, Y.; Li, P.; Zhang, Q.; Lei, J.; Zhang, Z.; Ding, X.; Zhou, H.; Zhang, W. Nanobody-based enzyme immunoassay for aflatoxin in agro-products with high tolerance to cosolvent methanol. *Anal. Chem.* 2014, *86*, 8873–8880. [CrossRef] [PubMed]
- Wang, Y.; Wang, H.; Li, P.; Zhang, Q.; Kim, H.J.; Gee, S.J.; Hammock, B.D. Phage-displayed peptide that mimics aflatoxins and its application in immunoassay. *J. Agric. Food Chem.* 2013, *61*, 2426–2433. [CrossRef] [PubMed]
- 57. Wang, J.; Bever, C.R.; Majkova, Z.; Dechant, J.E.; Yang, J.; Gee, S.J.; Xu, T.; Hammock, B.D. Heterologous antigen selection of camelid heavy chain single domain antibodies against tetrabromobisphenol A. *Anal. Chem.* **2014**, *86*, 8296–8302. [CrossRef] [PubMed]
- Pírez-Schirmer, M.; Rossotti, M.; Badagian, N.; Leizagoyen, C.; Brena, B.M.; Gonzalez-Sapienza, G. Comparison of three anti-hapten VHH selection strategies for the development of highly sensitive immunoassays for microcystins. *Anal. Chem.* 2017, *89*, 6800–6806. [CrossRef] [PubMed]
- 59. Doyle, P.J.; Arbabi-Ghahroudi, M.; Gaudette, N.; Furzer, G.; Savard, M.E.; Gleddie, S.; Mclean, M.D.; Mackenzie, C.R.; Hall, J.C. Cloning, expression, and characterization of a single-domain antibody fragment with affinity for 15-acetyl-deoxynivalenol. *Mol. Immunol.* **2008**, 45, 3703–3713. [CrossRef] [PubMed]
- 60. Anderson, G.P.; Goldman, E.R. TNT detection using llama antibodies and a two-step competitive fluid array immunoassay. *J. Immunol. Methods* **2008**, *339*, 47–54. [CrossRef] [PubMed]
- 61. Gonzalez-Sapienza, G.; Rossotti, M.A.; Rosa, S.T. Single-domain antibodies as versatile affinity reagents for analytical and diagnostic applications. *Front. Immunol.* **2017**, *8*, 977. [CrossRef] [PubMed]
- 62. Bever, C.R.; Majkova, Z.; Radhakrishnan, R.; Suni, I.; Mccoy, M.; Wang, Y.; Dechant, J.; Gee, S.; Hammock, B.D. Development and utilization of camelid VHH antibodies from alpaca for 2,2',4,4'-tetrabrominated diphenyl ether detection. *Anal. Chem.* **2014**, *86*, 7875–7882. [CrossRef] [PubMed]
- 63. Ahn, K.C.; Gee, S.J.; Tsai, H.J.; Bennett, D.; Nishioka, M.G.; Blum, A.; Fishman, E.; Hammock, B.D. Immunoassay for monitoring environmental and human exposure to the polybrominated diphenyl ether BDE-47. *Environ. Sci. Technol.* **2009**, *43*, 7784–7790. [CrossRef] [PubMed]
- 64. Xu, C.; Ou, J.; Cui, Y.; Wang, L.; Lv, C.; Liu, K.; Wang, B.; Xu, T.; Li, Q.; Liu, S. Development of a monoclonal antibody-based enzyme-linked immunosorbent assay for tetrabromobisphenol A. *Monoclon. Antib. Immunodiagn. Immunother.* **2013**, *32*, 113–118. [CrossRef] [PubMed]

- Xu, T.; Wang, J.; Liu, S.; Lv, C.; Shelver, W.L.; Li, Q.; Li, J. A highly sensitive and selective immunoassay for the detection of tetrabromobisphenol A in soil and sediment. *Anal. Chim. Acta* 2012, 751, 119–127. [CrossRef] [PubMed]
- Reddy, S.V.; Mayi, D.K.; Reddy, M.U.; Thirumala-Devi, K.; Reddy, D.V. Aflatoxins B1 in different grades of chillies (Capsicum annum L.) in India as determined by indirect competitive-ELISA. *Food Addit. Contam.* 2001, *18*, 553–558. [CrossRef] [PubMed]
- 67. Liu, X.; Tang, Z.; Duan, Z.; He, Z.; Mei, S.; Wang, X.; Gee, S.J.; Hammock, B.D.; Xu, Y. Nanobody-based enzyme immunoassay for ochratoxin A in cereal with high resistance to matrix interference. *Talanta* **2017**, *164*, 154–158. [CrossRef] [PubMed]
- 68. Li, X.; Li, P.; Zhang, Q.; Zhang, Z.; Li, R.; Zhang, W.; Ding, X.; Chen, X.; Tang, X. A sensitive immunoaffinity column-linked indirect competitive ELISA for ochratoxin A in cereal and oil products based on a new monoclonal antibody. *Food Anal. Methods* **2013**, *6*, 1433–1440. [CrossRef]
- 69. Thirumala-Devi, K.; Mayo, M.A.; Reddy, G.; Reddy, S.V.; Delfosse, P.; Reddy, D.V.R. Production of polyclonal antibodies against ochratoxin A and its detection in chilies by ELISA. *J. Agric. Food Chem.* **2000**, *48*, 5079–5082. [CrossRef] [PubMed]
- 70. Sinha, R.C.; Savard, M.E.; Lau, R. Production of monoclonal antibodies for the specific detection of deoxynivalenol and 15-acetyldeoxynivalenol by ELISA. *J. Agric. Food Chem.* **1995**, *43*, 1740–1744. [CrossRef]
- 71. Usleber, E.; Schneider, E.; Martlbauer, E.; Terplan, G. Two formats of enzyme immunoassay for 15-acetyldeoxynivalenol applied to wheat. *J. Agric. Food Chem.* **1993**, *41*, 2019–2023. [CrossRef]
- 72. Wesolowski, J.; Alzogaray, V.; Reyelt, J.; Unger, M.; Juarez, K.; Urrutia, M.; Cauerhff, A.; Danquah, W.; Rissiek, B.; Scheuplein, F. Single domain antibodies: Promising experimental and therapeutic tools in infection and immunity. *Med. Microbiol. Immunol.* **2009**, *198*, 157–174. [CrossRef] [PubMed]
- 73. Zhang, Q.; Wu, Y.; Wang, L.; Hu, B.; Li, P.; Liu, F. Effect of hapten structures on specific and sensitive enzyme-linked immunosorbent assays for N-methylcarbamate insecticide metolcarb. *Anal. Chim. Acta* 2008, 625, 87–94. [CrossRef] [PubMed]
- 74. Fodey, T.L.; Greer, N.M.; Crooks, S.R. Antibody production: Low dose immunogen vs. low incorporation hapten using salmeterol as a model. *Anal. Chim. Acta* 2009, *637*, 328–332. [CrossRef] [PubMed]
- 75. Zhang, D.; Li, P.; Yang, Y.; Zhang, Q.; Zhang, W.; Xiao, Z.; Ding, X. A high selective immunochromatographic assay for rapid detection of aflatoxin B1. *Talanta* **2011**, *85*, 736–742. [CrossRef] [PubMed]
- Li, P.; Wu, J.; Zhang, L.; Fan, Z.; Yu, T.; Jiang, F.; Tang, X.; Zhang, Z.; Zhang, W.; Zhang, Q. Doses of immunogen contribute to specificity spectrums of antibodies against aflatoxin. *Toxins* 2017, *9*, 172. [CrossRef] [PubMed]
- 77. van Houwelingen, A.; De Saeger, S.; Rusanova, T.; Waalwijk, C.; Beekwilder, J. Generation of recombinant alpaca VHH antibody fragments for the detection of the mycotoxin ochratoxin A. *World Mycotoxin J.* **2008**, *1*, 407–417. [CrossRef]
- Pérez, J.M.; Renisio, J.G.; Prompers, J.J.; van Platerink, C.J.; Cambillau, C.; Darbon, H.; Frenken, L.G. Thermal unfolding of a llama antibody fragment: A two-state reversible process. *Biochemistry* 2001, 40, 74–83. [CrossRef] [PubMed]
- 79. De Meyer, T.; Muyldermans, S.; Depicker, A. Nanobody-based products as research and diagnostic tools. *Trends Biotechnol.* **2014**, *32*, 263–270. [CrossRef] [PubMed]
- Muyldermans, S. Single domain camel antibodies: Current status. *Rev. Mol. Biotechnol.* 2001, 74, 277–302. [CrossRef]
- Liu, X.; Xu, Y.; Wan, D.; Xiong, Y.; He, Z.; Wang, X.; Gee, S.J.; Ryu, D.; Hammock, B.D. Development of a nanobody-alkaline phosphatase fusion protein and its application in a highly sensitive direct competitive fluorescence enzyme immunoassay for detection of ochratoxin A in cereal. *Anal. Chem.* 2015, *87*, 1387–1394. [CrossRef] [PubMed]
- 82. Liu, X.; Xu, Y.; Xiong, Y.; Tu, Z.; Li, Y.; He, Z.; Qiu, Y.; Fu, J.; Gee, S.J.; Hammock, B.D. VHH phage-based competitive real-time immuno-polymerase chain reaction for ultrasensitive detection of ochratoxin A in cereal. *Anal. Chem.* **2014**, *86*, 7471–7477. [CrossRef] [PubMed]
- Lei, J.; Li, P.; Zhang, Q.; Wang, Y.; Zhang, Z.; Ding, X.; Zhang, W. Anti-idiotypic nanobody-phage based real-time immuno-PCR for detection of hepatocarcinogen aflatoxin in grains and feedstuffs. *Anal. Chem.* 2014, *86*, 10841–10846. [CrossRef] [PubMed]

- Tang, X.; Li, P.; Zhang, Q.; Zhang, Z.; Zhang, W.; Jiang, J. Time-resolved fluorescence immunochromatographic assay developed using two idiotypic nanobodies for rapid, quantitative, and simultaneous detection of aflatoxin and zearalenone in maize and its products. *Anal. Chem.* 2017, *89*, 11520–11528. [CrossRef] [PubMed]
- Qiu, Y.L.; He, Q.H.; Xu, Y.; Wang, W.; Liu, Y.Y. Modification of a deoxynivalenol-antigen-mimicking nanobody to improve immunoassay sensitivity by site-saturation mutagenesis. *Anal. Bioanal. Chem.* 2016, 408, 895–903. [CrossRef] [PubMed]
- 86. Zarebski, L.M.; Urrutia, M.; Goldbaum, F.A. Llama single domain antibodies as a tool for molecular mimicry. *J. Mol. Biol.* **2005**, *349*, 814–824. [CrossRef] [PubMed]
- 87. Bazin, I.; Tria, S.A.; Hayat, A.; Marty, J.L. New biorecognition molecules in biosensors for the detection of toxins. *Biosens. Bioelectron.* **2017**, *87*, 285–298. [CrossRef] [PubMed]
- Xiong, Y.; Tu, Z.; Huang, X.; Xie, B.; Xiong, Y.H.; Xu, Y. Magnetic beads carrying poly(acrylic acid) brushes as "nanobody containers" for immunoaffinity purification of aflatoxin B1 from corn samples. *RSC Adv.* 2015, *5*, 77380–77387. [CrossRef]
- Li, X.; Li, P.; Zhang, Q.; Li, Y.; Zhang, W.; Ding, X. Molecular characterization of monoclonal antibodies against aflatoxins: A possible explanation for the highest sensitivity. *Anal. Chem.* 2012, *84*, 5229–5235. [CrossRef] [PubMed]
- Spinelli, S.; Frenken, L.G.; Hermans, P.; Verrips, T.; Brown, K.; Tegoni, M.; Cambillau, C. Camelid heavy-chain variable domains provide efficient combining sites to haptens. *Biochemistry* 2000, *39*, 1217–1222. [CrossRef] [PubMed]
- 91. Spinelli, S.; Tegoni, M.; Frenken, L.; Van, V.C.; Cambillau, C. Lateral recognition of a dye hapten by a llama VHH domain. *J. Mol. Biol.* **2001**, *311*, 123–129. [CrossRef] [PubMed]
- 92. Fanning, S.W.; Horn, J.R. An anti-hapten camelid antibody reveals a cryptic binding site with significant energetic contributions from a nonhypervariable loop. *Protein Sci.* 2011, 20, 1196–1207. [CrossRef] [PubMed]
- 93. Singh, A.; Pasha, S.K.; Manickam, P.; Bhansali, S. Single-domain antibody thermally stable electrochemical immunosensor. *Biosens. Bioelectron.* **2016**, *83*, 162–168. [CrossRef] [PubMed]
- 94. Xu, L.; Zhang, Z.W.; Zhang, Q.; Li, P.W. Mycotoxin determination in foods using advanced sensors based on antibodies or aptamers. *Toxins* **2016**, *8*, 239. [CrossRef] [PubMed]
- Chen, A.; Wang, R.; Bever, C.R.S.; Xing, S.; Hammock, B.D.; Pan, T. Smartphone-interfaced lab-on-a-chip devices for field-deployable enzyme-linked immunosorbent assay. *Biomicrofluidics* 2014, *8*, 3187–3194. [CrossRef] [PubMed]



© 2018 by the authors. Licensee MDPI, Basel, Switzerland. This article is an open access article distributed under the terms and conditions of the Creative Commons Attribution (CC BY) license (http://creativecommons.org/licenses/by/4.0/).