

Supplementary material: Headspace solid-phase microextraction gas chromatography-mass spectrometry analysis of scent profiles from human skin

Emer Duffy, Gwendoline Albero and Aoife Morrin




Figure S1. Distribution of compound classes in human skin volatile samples. Compounds were classified according to chemical functionality based on chromatographic peak areas as determined by HS-SPME GC-MS. (F – female, M – male)




Figure S2. Box plot comparing intra- and inter-participant variance for hexanal.




Figure S3. Box plot comparing intra- and inter-participant variance for octanal.




Figure S4. Box plot comparing intra- and inter-participant variance for nonanal.




Figure S5. Box plot comparing intra- and inter-participant variance for *n*-decanoic acid.




Figure S6. Box plot comparing intra- and inter-participant variance for tetradecanoic acid.




Figure S7. Box plot comparing intra- and inter-participant variance for *n*-hexadecanoic acid.




Figure S8. Box plot comparing intra- and inter-participant variance for isopropyl palmitate.




Figure S9. Box plot comparing intra- and inter-participant variance for squalane.




Figure S10. Box plot comparing intra- and inter-participant variance for 6-methyl-5-hepten-2-one.




Figure S11. Box plot comparing intra- and inter-participant variance for geranyl acetone.

Table S1. Measured tissue dielectric constants (TDC) ($\pm 95\%$ confidence intervals) for participants' skin (F – female, M – male, $n = 3$).

Participant	TDC
F1	24.0 ± 1.0
F2	26.4 ± 0.3
F3	25.4 ± 0.7
F4	29.0 ± 1.2
M1	34.5 ± 1.4
M2	42.7 ± 1.5
M3	35.9 ± 1.6
M4	39.0 ± 1.2




Figure S12. Participant scent profiles after addition of fragrance to skin. Percentage contribution of olfactory family is based on peak areas obtained by headspace solid-phase microextraction sampling with gas chromatography-mass spectrometry analysis. (Y-axis labels show participant and time of sampling *e.g.* F1_0 represents Female 1 sample collected at 0 hr).