

Article



Determination of Parabens and Phenolic Compounds in Dairy Products through the Use of a Two-Step Continuous SPE System Including an Enhanced Matrix Removal Sorbent in Combination with UHPLC—MS/MS

Laura Palacios Colón, Andrés J. Rascón 🗅 and Evaristo Ballesteros *🕩

Department of Physical and Analytical Chemistry, E.P.S of Linares, University of Jaén, 23700 Linares, Jaén, Spain; lpcolon@ujaen.es (L.P.C.); arascon@ujaen.es (A.J.R.) * Correspondence: eballes@ujaen.es

Abstract: Dairy products can be contaminated by parabens and phenolic compounds from a vast variety of sources, such as packaging and manufacturing processes, or livestock through feed and environmental water. A two-step continuous solid-phase extraction (SPE) and purification methodology was developed here for the determination of both types of compounds. In the first step, a sample extract is passed in sequence through an EMR-lipid sorbent and an Oasis PRiME HBL sorbent to remove fat and preconcentrate the analytes for subsequent detection and quantification by UHPLC–MS/MS. This method enabled the determination of 28 parabens and phenolic contaminant with excellent recovery (91–105%) thanks to the SPE sorbent combination used. The proposed method was validated through the determination of the target compounds, and was found to provide low detection limits (1–20 ng/kg) with only slight matrix effects (0–10%). It was used to analyse 32 different samples of dairy products with different packaging materials. Bisphenol A and bisphenol Z were the two phenolic compounds quantified in the largest number of samples, at concentrations over the range of 24–580 ng/kg, which did not exceed the limit set by European regulations. On the other hand, ethylparaben was the paraben found at the highest levels (33–470 ng/kg).

Keywords: phenolic compounds; parabens; dairy products; continuous solid-phase extraction; EMR-lipid sorbent; UHPLC-MS/MS

1. Introduction

Dairy products are highly complex foods rich in micro- and macro-nutrients that are massively consumed in Spain. Based on the latest official food consumption report, milk is consumed at a rate of 71.27 L/person each year in the country [1]. Milk and its derivatives are exposed to contamination from the environment, livestock feed, and packaging materials. Environmental pollution is a source of major concern as it jeopardizes the sustainability of the planet. Industrialization, the rapid growth of the human population, and anthropogenic activities generate contamination that pollutes the ecosystem [2].

Parabens, phenolic compounds, and various other pollutants can reach milk through cattle (cows, goats, and sheep) being directly exposed to environmental or feeding contaminants [3]. The long-term consumption of foods of an animal origin containing residual contaminants can have toxic effects on the human body through bioaccumulation [4]. Some compound families have been deemed endocrine disruptor chemicals (EDCs) because they can interfere with biosynthetic processes and hormonal action. Most EDCs' phenolic compounds and the bisphenols included come from anthropogenic sources [5,6].

Parabens are alkyl esters derived from *p*-hydroxybenzoic acid, used by the cosmetic, pharmaceutical, food, and beverage industries for their antimicrobial and preservative properties [6–8]; however, phenolic compounds are highly toxic and persistent in the



Citation: Palacios Colón, L.; Rascón, A.J.; Ballesteros, E. Determination of Parabens and Phenolic Compounds in Dairy Products through the Use of a Two-Step Continuous SPE System Including an Enhanced Matrix Removal Sorbent in Combination with UHPLC–MS/MS. *Foods* **2023**, *12*, 2909. https://doi.org/10.3390/ foods12152909

Academic Editor: Salvador Maestre Pérez

Received: 30 June 2023 Revised: 19 July 2023 Accepted: 28 July 2023 Published: 31 July 2023



Copyright: © 2023 by the authors. Licensee MDPI, Basel, Switzerland. This article is an open access article distributed under the terms and conditions of the Creative Commons Attribution (CC BY) license (https:// creativecommons.org/licenses/by/ 4.0/). environment, so they pose some risks for living organisms. Phenolics are essential in the production of pesticides, pharmaceuticals, packaging materials, household preservatives, and as intermediates in industry [9]. For example, triclosan is usually added to body care products such as liquid and solid soap, toothpaste, and mouthwash. Its antimicrobial properties also make it useful for food contact materials, where it can account for as much as 1% (w/w) of their contents, so it has been banned for use in medical devices, for example [10]. Bisphenol A (BPA) is a synthetic phenol used to obtain polycarbonate plastics and epoxy resins for the production of a number of consumer products [8,11,12]. BPA is suspected to cause health issues such as cancer and reproductive disorders, and also to be an obesity promoter [9,13]. On these grounds, the European Food Safety Authority (EFSA) recently downgraded its tolerable daily intake (TDI) 100,000 times (from 4000 ng/kg to 0.04 ng/kg) [14]. In contrast, the European Commission set the maximum specific migration limit (SML) for BPA at 0.05 mg/kg in 2018 to update the previous 2011 level [15]. BPA has been already banned in food contact materials intended for children, which has forced manufacturers to use alternative compounds such as bisphenol B, F, or S. As a result, these BPA analogues, which are similarly toxic to the parent compound or even more so, are now ubiquitous in environmental and biological samples, and the exposure of humans (through the consumption of dairy products, for example), other living organisms, and the environment has increased with time [16]. Because of its similarity in structure, toxicity, and metabolism, BPF has been deemed an unsafe alternative to BPA [17]. So far, the European Union has set the SML for only one BPA analogue, bisphenol S, at 0.05 mg/kg [18].

The US Environmental Protection Agency (EPA) know about the potential threat of phenol pollutants and has included the following to its Priority Pollutant List: 4-chloro-3-methylphenol, 2,4,6-trichlorophenol, 2,4-dichlorophenol, pentachlorophenol, 2,4-dinitrophenol, 2-methyl-4,6-dinitrophenol, 4-nitrophenol, phenol, 2-nitrophenol, and 2,4-dimethylphenol [19]. In addition, it has classified pentachlorophenol as a group B2, probable human carcinogen [13]. On the other hand, European Commission (EC) Regulation No. 10/2011 has set the SML for phenolics such as 4-*tert*-butylphenol, phenol and 2,6-dimethylphenol at 0.05 mg/kg for specific migration from food contact materials [18]. Nonylphenol and 4-octylphenol can disrupt reproductive and immune functions, and have adverse health effects as a result [12]. The former has been included as a priority contaminant on EPA's Candidate Contaminant List (CCL4) [20].

Parabens and phenolic compounds in dairy products are most often determined by gas or liquid chromatography. Gas chromatography (GC) is an effective choice for determining EDCs by virtue of its good separation efficiency and high throughput. For example, González et al. [21] and Palacios et al. [22] determined ECDs in dairy products by using GC coupled to a mass spectrometer (GC–MS), whereas Zhang et al. determined the phenolics in milk by GC with flame ionization detection (GC–FID) [23]. These compounds are also frequently determined by high-performance liquid chromatography (HPLC). The main advantage of HPLC over GC for phenolics is that the analytes need not be derivatized. Santonicola et al. [17], Tuzimski and Szubartowski [16], and Xiong et al. [24] used HPLC with fluorescence detection (HPLC–FLD) to determine the bisphenols in dairy products. However, as can be seen in Table 1, liquid chromatography coupled to tandem mass spectrometry (HPLC–MS/MS) is the most frequently employed analytical technique to monitor phenol compounds, bisphenols, and parabens, a result of its affording detection at trace levels in dairy products [7,25–28].

Sample pre-treatment is a fundamental step in the determination of the previous types of compounds, especially when analytes are present at low concentrations in complicated matrices such as dairy products. Solid-phase extraction (SPE) provides a highly efficient, simple, and economical method to enrich samples with the target analytes for their subsequent determination [29]. As can be seen in Table 1, SPE has been used with various sorbent materials to extract parabens and phenolic pollutants from milk and dairy products [17,28,30]. In addition, a new generation of easily used, robust sample preparation sorbents with virtually no matrix effects, including Oasis PRiME HLB, has recently been

developed [31]. This is a reversed-phase SPE sorbent that is packed in a column where the sample is directly passed through after extraction, saving time, and it is also more environmentally friendly than the conventional procedure and can remove the majority of matrix interferents such as proteins, salts, and phospholipids [4]. Enhanced matrix removal-lipid (EMR-lipid), which was introduced by Agilent Technology in 2015, is another recently developed sorbent material especially designed for lipid removal. In fact, EMR-lipid selectively clear major lipids from samples avoiding undesired analyte retention [16,27,32,33].

As can be seen in Table 1, there are other techniques in the literature for the treatment of milk and dairy product samples, such as QuEChERS [24,25,27], matrix solid-phase dispersion [7], dispersive magnetic micro-solid phase extraction [34], electrochemical assistance solid-phase microextraction [23], and ultrasound-assisted liquid—liquid microextraction [35].

In previous work, our group determined EDCs using GC-MS [22], which required derivatizing the analytes to boost analytical signals and avoid false positives. In this work, we determined seven parabens and twenty-one phenolic compounds classified as EDCs in dairy products by using ultra-high-performance liquid chromatography-tandem mass spectrometry (UHPLC–MS/MS) to expedite sample treatment and avoid the need for organic solvents. This technique is one of the best choices for quantitative analysis thanks to its high sensitivity and reproducibility, confirmed by the two most abundant product ions. EMR-lipid was used in combination with Oasis PRiME HLB in a closed semi-automated SPE system to remove fat and protein from the sample matrix and increase the recovery of the target analytes as a result. The main goal was to simplify and unify the extraction procedure for the target compounds while improving the efficiency of the process and reducing the matrix effects. The resulting SPE-UHPLC-MS/MS method will be evaluated in terms of sensitivity, selectivity, precision, and recovery, after which it will be used to determine the target analytes in dairy products packed in different kinds of food contact materials (namely, high density polyethylene (HDPE), polystyrene (PS), polyethylene terephthalate (PET), multilayer packaging, and glass) to check whether any of them were responsible for the presence of the analytes in any sample.

Analysias a	Samulas	Samula Treatment b	Tachniqua	Analytical Figures of Marit d	Concentrations Found in Peol Formulas	Poforoncos
Analytes ~	Samples	Sample Treatment	lecnnique ·	I OO 0 91–4 4 ng/kg	Concentrations Found in Real Samples	Kererences
BPA, parabens, and other plasticisers	Milk	MSPD	HPLC-MS/MS	RSD < 14%	BPA 1.8–59 ng/kg	[7]
BPF, BPE, BPB, BPA, and diglycidyl ether	Breast milk	SPE, QuEChERS	HPLC-FLD	LOD 57,000–78,000 ng/L RSD 1.0–14% R 57–88%	n.d.	[16]
BPA and BPF	Milk	SPE	HPLC-FLD	LOD 30 ng/L RSD 1.0-14.8% R 98-107%	BPF 1020–2690 ng/L	[17]
BPA, BPB, BPF, BPZ, BPS, and other bisphenol analogues	Yogurt	QuEChERS	CG-MS	LOQ 170 ng/kg	n.d.	[21]
Phenol, alkyphenols, chlorophenols, and bisphenols	Dairy products	SPE	CG/MS	LOD 6–63 ng /L RSD 2.4–11% R 85–108%	BPA 30–1400 ng/kg BPZ 96–1100 ng/kg BPF 270–950 ng/kg Nonylphenol 56–390 ng/kg 4-tert-Butylphenol 310–2100 ng/kg 3.4-Dimethylphenol 130–2300 ng/kg 4-Pentylphenol 190–990 ng/kg 2-Phenylphenol 320 ng/kg	[22]
4-Chlorophenol, 4- <i>tert</i> -butylphenol, 4- <i>tert</i> -octylphenol, and other phenols	Milk	EA-SPME	GC-FID	LOD 1–30 ng/L RSD 3.8–12.4% R 88–119%	4-tert-Butylphenol 850 ng/L 4-Chlorophenol n.d. 4-tert-Octylphenol n.d.	[23]
BPA, BPB, and other bisphenol analogues	Milk	QuEChERS	HPLC-FLD	LOD 1000–3100 ng/L RSD 9.1–16% R 86–99%	BPA 1374 ng/kg	[24]
BPA, BPF, BPS, and parabens	Breast milk	QuEChERS	HPLC-MS/MS	LOQ 100–250 ng/L RSD 5–16% R 83–115%	BPA 130–1 620 ng/L BPF 130–320 ng/L BPS 370 ng/L Parabens 130–4050 ng/L	[25]
BPA, nonylphenol, and 4-ter-octylphenol	Milk	MSPD	HPLC-MS/MS	LOD 50–100 ng/kg RSD 1.4–7.8% R 84–103%	BPA 490 ng/kg Nonylphenol 100 ng/kg 4- <i>ter-</i> Octylphenol 4240–17,600 ng/kg	[26]
BPA, BPB, BPF, BPS, and other bisphenol analogues	Breast milk and sweetened condensed milk	QuEChERS	HPLC-DAD/ HPLC-MS/MS	LOQ 100–250 ng/L RSD 8–17% R 35–102%	BPA 230–690 ng/L BPS 290–680 ng/L BPF 220–290 ng/L	[27]
Parabens, TCS, BPA, nonylphenol, 4- <i>tert</i> -octylphenol, and other phenols	Breast milk	SPE	CG/MS	LOD 1–9 ng /L RSD 4.4–7.0% R 86–104%	Phenols 550–5600 ng/L BPA 1400–2900 ng/L Parabens 15–8100 ng/L	[30]
BPA	Milk	DMSPE	HPLC-UV	LOD 3 050 ng/L RSD 9.1–16% R 86–99%	n.d. ^e	[34]

Table 1. Selected references for the determination of parabens and phenolic compounds in dairy products.

Table 1. Cont.

Analytes ^a	Samples	Sample Treatment ^b	Technique ^c	Analytical Figures of Merit ^d	Concentrations Found in Real Samples	References
BPA, BPF, BPAF, and 4-chlorophenol	Milk	UA-DLLME	HPLC-UV	LOD 250–1000 ng/L RSD 0.67–13.7% R 83–119%	n.d.	[35]
BPA and BPS	Milk	LLE	PSI-MS	LOD 1500–4800 ng /L RSD 5.9–18.9% R 74–112%	BPA 77,600–150,800 ng/L BPS 60,000 ng/L	[36]
Nonylphenol, 4-tert-octylphenol and BPA	Dairy products	UA-DLLME	HPLC-FLD	LOD 10–40 ng/kg RSD 2.8–8.8% R 83–112%	BPA 800–128,700 ng/kg Nonylphenol 3500–36,700 ng/kg 4- <i>tert</i> -Octylphenol 300–29,500 ng/kg	[37]
BPA, BPB, BPF, BPZ, 4- <i>tert</i> -octylphenol, and nonylphenol	Milk	QuEChERS	UHPLC-MS/MS HPLC-FLD	LOD 200–1500 ng/L	BPA 500–5600 ng/L BPF 500–8700 ng/L	[38]
BPA, BPS, BPAF, BPF and parabens	Breast milk	QuEChERS	HPLC-MS/MS	LOD 10–200 ng/L RSD 1.0–7.9% R 77–100%	BPA 990–1910 ng/L BPS 60–200 ng/L BPAF 60–70 ng/L BPF n.d. Parabens 30–1520 ng/L	[39]
Parabens, phenol, alkylphenols, chlorophenols and bisphenols	Dairy products	Continuous SPE	HPLC-MS/MS	LOD 1–20 ng/kg RSD 2.0–11.4% R 91–105%	BPA 33–580 ng/kg BPZ 24–57 ng/kg Nonylphenol 34–210 ng/kg 4- <i>tert</i> -Butylphenol 21–220 ng/kg 4-Chlorophenol 35–160 ng/kg Pentachlorophenol 35–76 ng/kg Parabens 9–470 ng/kg	This work

^a BPA, bisphenol A; BPAF, bisphenol AF; BPB, bisphenol B; BPE, bisphenol E; BPF, bisphenol F; BPS, bisphenol S; BPZ, bisphenol Z. ^b DMSPE, dispersive magnetic solid-phase extraction; EA–SPME, electrochemical assistance solid-phase microextraction; MSPD, matrix solid-phase dispersive liquid–liquid microextraction; VASUSME, vortex-assisted supramolecular solvent microextraction; SPME, solid-phase micro-extraction; UA–DLLME, ultrasound–assisted dispersive liquid–liquid microextraction; VASUSME, vortex-assisted supramolecular solvent microextraction; HPLC–DAD, high-performance liquid chromatography–photodiode array detection; HPLC–FLD, high-performance liquid chromatography with fluorescence detection; HPLC–MS/MS, high-performance liquid chromatography coupled to triple-quadrupole mass spectrometry; HPLC–UV, high-performance liquid chromatography with ultraviolet detection; GC–FID, gas chromatography with flame–ionization detection; GC–MS, gas chromatography coupled to single-quadrupole mass spectrometry; MSPD, matrix solid-phase dispersion; PSI-MS, paper spray ionization–mass spectrometry. ^d LOD, limit of detection; LOQ, limit of quantification; R, recovery; RSD, relative standard deviation. ^e n.d., not detected.

2. Materials and Methods

2.1. Reagents and Chemicals

HPLC-grade water, methanol (MeOH), acetone, acetonitrile (ACN), and *n*-hexane were supplied by Sigma-Aldrich (St. Louis, MO, USA). Formic acid (FA, 98% purity) was purchased from Fluka (St. Louis, MO, USA), and MgSO₄ and NaCl were both obtained in the highest available purity from PanReac AppliChem (Barcelona, Spain).

Analytical standards of parabens (metylparaben, ethylparaben, propylparaben, isopropylparaben, butylparaben, isobutylparaben, and benzylparaben), phenolic compounds (3,4-dimethylphenol, 2,5-dimethylphenol, 4-*tert*-butylphenol, 2-*tert*-butyl-4-methylphenol, 4pentylphenol, 2-phenylphenol, 4-phenylphenol, 4-hexylphenol, 4-heptylphenol, 4-chlorophenol, 4-chloro-3-methylphenol, pentachlorophenol, triclosan, 4-*tert*-octylphenol, nonylphenol, BPA, bisphenol F, bisphenol B, bisphenol S, and bisphenol Z), and triphenylphosphate (internal standard (IS)) were supplied by Sigma-Aldrich (St. Louis, MO, USA) in the highest available purity.

Centrifuge tubes were obtained from J.D. Catalan S.L. (Madrid, Spain). Regarding the sorbents, Oasis HLB and Oasis PRiME HLB were acquired from Waters Corporation (Milford, MA, USA), and Captiva EMR-Lipid was acquired from Agilent Technologies (Santa Clara, CA, USA).

All of the analytical standards were made using 5 g/L in MeOH. The standard solutions were stored in glass-stoppered bottles at 4 °C in the dark. Mixed standards with all of the analytes were prepared at 1 mg/L in MeOH on a daily basis and diluted when needed. The SPE eluent was 30:70 (v/v) ACN:H₂O containing of 400 µg/L of IS and was prepared daily.

2.2. Equipment

For separation, a Dionex Ultimate 3000 UHPLC chromatograph (Thermo Fisher Scientific, Waltham, MA, USA) was used equipped with Zorbax Rapid Resolution High Definition (RRHD) Eclipse Plus C18 column (2.1 mm \times 100 mm, 1.8 µm particle size) from Agilent Technologies (Santa Clara, MA, USA). The mobile phases were ultrapure water (Solvent A) and acetonitrile (Solvent B), both containing 0.1% of formic acid; the flow was set at 0.4 mL/min. The column oven was set at 40 °C and the injection volume was set to 10 µL. The elution programme was as follows: 10% B (0–2 min), 10–30% B (2–18 min), 30–40% B (18–24 min), 40–70% B (24–32 min), and 70–10% B (32–35 min). The UHPLC equipment was coupled to a TSQ Endura triple quadrupole (QqQ) from Thermo Scientific (USA) equipped with a heated electrospray ionization probe (HESI) set in the positive ion mode. The spray voltage was 3500 V; the sheath and aux gas settings were 45 and 13 arbitrary units, respectively; the ion transfer tube and vaporizer temperature were 400 and 350 °C, respectively; and the collision gas (CID) pressure was 2 kPa. The individual analytes were infused and characterized in the mass spectrometer triple quadrupole equipment in multiple reaction monitoring mode (MRM) (Table 2).

	Compounds	t _R ^a	Precursor (m/z)	Quantification Peak (<i>m</i> / <i>z</i>)	Collision Energy (V)	Confirmation Peak (<i>m</i> / <i>z</i>)	Collision Energy (V)	RF Lens (V)	LOD ^b (ng/kg)	Linearity (ng/kg)	r ^{2 c}
	Methylparaben	4.63	152.2	93.0	21.0	137.0	14.1	56.0	1.0	3.3-20,000	0.997
sua	Ethylparaben	6.57	166.2	93.0	22.6	138.0	14.0	50.0	1.0	3.4-20,000	0.996
	Isopropylparaben	9.48	179.3	93.0	21.6	137.0	14.3	49.6	1.0	3.3-20,000	0.994
abe	Butylparaben	9.66	197.0	129.1	11.7	178.1	12.5	51.8	1.1	3.6-20,000	0.996
ar	Propylparaben	10.11	180.2	93.0	23.3	92.0	23.8	53.8	1.1	3.6-20,000	0.995
H	Isobutylparaben	15.29	193.2	92.0	24.8	137.1	15.7	51.3	10	34-20,000	0.997
	Benzylparaben	16.42	227.1	92.0	24.1	136.1	14.9	50.1	1.0	3.3–20,000	0.998
	Pentylphenol	0.74	165.0	96.9	10.3	78.9	26.7	37.4	1.1	3.5-20,000	0.997
	4-Chloro-3-methylphenol	0.97	141.0	35.0	19.1	105.0	16.7	62.4	10	34-20,000	0.994
	4-tert-Octylphenol	1.15	209.1	96.9	15.0	78.9	34.1	65.8	10	35-20,000	0.994
	Phenol	1.27	91.1	45.0	10.3	89.4	55.0	206.5	10	34-20,000	0.999
	4-tert-Butylphenol	1.48	147.1	61.9	11.4	46.0	45.2	35.7	1.1	3.5-20,000	0.995
	Nonylphenol	1.59	219.1	172.9	10.3	171.0	19.9	80.9	1.1	3.6-20,000	0.998
	2-tert-Butyl-4-methylphenol	1.69	147.1	61.9	13.9	89.0	12.8	30.0	1.0	3.3-20,000	0.997
	Bisphenol F	1.82	198.0	130.1	10.3	171.1	10.8	30.0	10	34-20,000	0.994
	3.4-Dimethylphenol	3.44	121.0	77.0	13.1	91.9	23.5	52.8	20	65-20,000	0.995
ols	2.5-Dimethylphenol	3.66	121.0	76.9	10.9	91.9	25.8	49.8	10	33-20,000	0.997
enc	4-Phenylphenol	4.19	169.0	19.2	32.6	69.0	54.1	38.7	1.0	3.3-20,000	0.995
Ρh	2-Phenylphenol	4.36	169.1	19.3	41.9	150.9	13.8	48.3	1.1	3.6-20,000	0.997
	Bisphenol S	4.51	249.1	108.0	28.0	92.0	35.3	109.1	1.0	3.4-20,000	0.997
	4-Chlorophenol	7.11	127.0	35.1	19.3	91.0	17.6	60.4	7.0	23-20,000	0.994
	4-Hexylphenol	16.74	177.1	148.0	16.0	149.1	16.8	158.8	10	34-20,000	0.995
	Bisphenol A	16.89	221.1	205.2	30.2	164.0	25.2	115.8	5.1	17-20,000	0.996
	Bisphenol B	19.08	241.1	225.0	10.2	223.0	10.2	68.6	5.0	16-20,000	0.994
	4-Heptylphenol	19.54	195.3	179.1	19.9	167.0	16.0	78.7	7.0	23-20,000	0.996
	Bisphenol Z	20.48	267.1	173.1	27.0	223.0	32.6	100.7	5.0	17-20,000	0.994
	Pentachlorophenol	25.8	264.8	35.0	55.0	97.0	27.4	125.4	10	34-20,000	0.994
	Triclosan	25.36	289.0	35.04	10.3	37.1	10.2	45.1	1.0	3.3–20,000	0.996
	Triphenylphosphate (IS)	27.54	337.3	96.9	32.7	79.9	46.5	122.9			

Table 2. UHPLC-MS/MS parameters for the determination of the studied da	airy product samples.
---	-----------------------

t_R^a: retention time; ^b LOD: limit of detection; ^c r²: correlation coefficient.

Other equipment used included the continuous solid-phase extraction system, composed of a Gilson peristaltic pump (Villiers-le-Bel, France), three injection valves of model 5041 from Rheodyne (Cotati, CA, USA), poly (vinyl chloride) tubing, PTFE columns that were custom-packed with the sorbent materials (Oasis PRiME HLB and EMR-lipid), and a Centrofriger BL-II centrifuge from JP Selecta (Barcelona, Spain). The Oasis PRiME HLB sorbent column was conditioned by passing 1 mL of methanol, 1 mL of acetonitrile, and 2 mL of Milli-Q water in this sequence. The EMR-lipid sorbent column was conditioned by passing 2 mL of Milli-Q water.

2.3. Sampling

Thirty-one samples of dairy products including fifteen samples of milk (skimmed, semi-skimmed, and whole milk from cow; fresh cow milk; semi-skimmed goat milk; and semi-skimmed sheep milk) and sixteen dairy products samples (yogurts, custards, milk-shakes, cheeses, creams, butters, and margarines) were analysed. The samples represent a wide range of the most consumed dairy products in Spain and were obtained from different commercial brands in containers made of various materials such as polyethylene terephthalate (PET), high-density polyethylene (HDPE), polystyrene (PS), multilayer material (a mixture of cardboard, polyethylene (PE) and aluminium), and glass. All of the samples were subjected to pasteurization, UHT, or sterilization for preservation. In addition, all were bought from Spanish supermarkets, kept under the conditions recommended on their packaging, and stored in the dark at 4 °C until analysis.

2.4. Sample Treatment

The sample pre-treatment of high-throughput screening methods usually involves retaining and extracting the target compounds as far as possible while removing impurities such as protein and fat to the greatest extent [4]. In this work, we developed an extraction procedure involving liquid—liquid extraction followed by clean-up and preconcentration of the extract with two solid-phase sorbents (EMR-lipid and Oasis PRiME HLB) in a serially connected semiautomated SPE system. The analytes were determined by UHPLC—MS/MS.

First, the samples were homogenized in 50 mL propylene centrifuge tubes. Then, 1 g of each sample (milk, milkshakes, custards, yogurts, cream, cheese, margarines, and butters) was transferred to another tube and supplied with 8 mL of ACN containing 0.5% FA (v/v), as previously described by our group [22]. The tubes were then centrifuged at 5000 rpm (2150 g) at 4 °C for 10 min, and the resulting supernatant was collected and supplied with 1.6 g MgSO₄ and 0.4 g NaCl, vortexed for 2 min, and centrifuged at 5000 rpm (2150 g) at 4 °C for 2 min.

The final extract was diluted to 50 mL with ultrapure water and its pH adjusted to 4 for passing through the closed, semi-automated SPE system packed with two serially columns of sorbent material. The first column contained 50 mg of EMR-lipid and the second 70 mg of Oasis PRiME HLB. The extract was introduced into the system at a flow rate of 5.0 mL/min to clean it up and allow the analytes to be retained on the sorbents, respectively, while the matrix and interferents were sent to waste (Figure 1). The Oasis PRiME HLB sorbent was then dried by an air stream circulating at a flow of 5 mL/min in both directions for 3 min and the analytes were then eluted with 600 μ L of a 30:70 (v/v) ACN:H₂O mixture containing a 400 μ g/L concentration of IS. The vials were hermetically sealed and a volume of 10 μ L was injected into the UHPLC–MS/MS system for analysis. Figure S1 (Supplementary Material) details the operation of the continuous system for the SPE.

cheese, and butter



Figure 1. Procedure for determining parabens and phenolic compounds in dairy products using UHPLC–MS/MS.

3. Results and Discussion

3.1. Sample Treatment Optimization

The good performance of a method for determining analytes at trace levels is largely dependent on how the samples are pre-treated. The matrices of dairy products are very complex due to their high fat (0.3–82%) and protein content (0.5–24%), thus requiring very efficient treatment procedures [40].

We initially used a liquid–liquid extraction (LLE) procedure previously optimized by our group [22]. For this purpose, an amount of 1 g of sample was spiked with a 500 ng/kg concentration of each analyte and subjected to LLE with 8 mL of acetonitrile containing 0.5% FA (v/v). SPE, which is one of the most commonly used techniques for the clean-up of extracts and for the preconcentration of analytes [17,41]. Because the analytes were highly soluble in fats, their presence in the extracts could detract from recoveries as the analytes might be co-discharged to waste with the fat fraction. The efficiency with which fat and other interferents were removed from the samples was assessed by processing the supernatant from the first step in the following ways: (1) direct analysis of the extract, (2) passage through an SPE system with the EMR-lipid sorbent, and (3) passage through an SPE system using a combination of EMR-Lipid and Oasis PRiME HBL. The efficiency of the sorbent was evaluated in terms of recovery. As can be seen in Figure 2, recoveries without EMR-lipid ranged from 20 to 120% for milk and 60 to 200% for butter. Passing the sample through the column containing EMR-lipid improved the recoveries somewhat (62–117% from milk and 42–163% from butter), although inadequately, and required introducing a third cleaning/preconcentration step involving passing the extract through a second column packed with Oasis PRiME HLB. Using the two sorbents substantially improved recoveries, which ranged from 90 to 105% for both milk and butter. Similar results were obtained with the other types of dairy products. Therefore, all further tests were conducted using a column packed with EMR-Lipid and then one packed with Oasis PRiME HBL.

The optimum amount of EMR-Lipid sorbent to be used to remove fat from the LLE extract was established by passing it through columns containing variable amounts of sorbent (20–80 mg). Amounts exceeding 40 mg were found to provide near-quantitative analyte recoveries. This led us to choose 50 mg of EMR-Lipid for further testing. The optimum amount of the second sorbent, Oasis PRiME HBL, was established similarly. A number of calibration graphs were run with different sorbent columns packed with 20–120 mg of sorbent to pass 100 mL of standard aqueous solutions with variable concentrations of each target compound over the range of 50–500 ng/L. The peak areas of the analytes increased with the increased amount of sorbent up to 60 mg, but decreased above 80 mg owing to the need for a higher volume of eluent to ensure complete elution of the analytes. Based on these results, columns packed with 70 mg of Oasis PRiME HBL were selected for further tests.







Figure 2. Recovery of parabens and phenolic compounds from milk (**A**) and butter (**B**) using different sorbents: no EMR-lipid, EMR-lipid, Oasis PRiME HLB, and EMR-lipid + Oasis PRiME HLB. (1) Phenol, (2) 2,5-dimethylphenol, (3) 4-chlorophenol, (4) 3,4-dimethylphenol, (5) 4-chloro-3-methylphenol, (6) 4-tert-butylphenol, (7) 2-tert-butyl-4-methylphenol, (8) metylparaben, (9) 4pentylphenol, (10) ethylparaben, (11) isopropylparaben, (12) 2-phenylphenol, (13) 4-hexylphenol, (14) 4-tert-octylphenol, (15) propylparaben, (16) isobutylparaben, (17) 4-heptylphenol, (18) nonylphenol, (19) butylparaben, (20) 4-phenylphenol, (21) pentachlorophenol, (22) triclosan, (23) bisphenol F, (24) benzylparaben, (25) bisphenol A, (26) bisphenol B (27) bisphenol Z, and (28) bisphenol S.

On the other hand, several solvents or a mixture of solvents (MeOH, EtOH, ACN, acetone, MeOH:H₂O, and ACN:H₂O) were tested as eluents of the analytes adsorbed in the continuous SPE system. For this purpose, 1 g of milk, milkshake, custard, heavy cream, yogurt, cheese, butter, or margarine was spiked with a 500 ng/kg concentration of each analyte and subjected to the treatment described in Section 2.4. The best results were obtained using the ACN:H₂O mixture, at a ratio of 30:70 (v/v), ensuring near-quantitative elution against only 10–40% with the other choices, also the effect of the volume used on the elution was examined in the range of 50–1000 µL. It was observed that for volumes

above 550 μ L, the elution of the analytes retained on the Oasis PRiME HBL column was quantitative. An eluent volume of 600 μ L was thus selected for further testing. It was also studied whether the sample flow rate affected the retention of analytes in the sorbent column. It could be seen that the sample flow rate that passed through the column had no effect on the sorption efficiency in the studied range (0.5–6.0 mL/min). Thus, 5.0 mL/min was selected as the sample and air flow rate. Finally, the sample breakdown volume was tested by using a determined amount of each compound (20 ng) in different volumes from 10 to 500 mL. Volumes up to 200 mL resulted in near-quantitative sorption with both Oasis PRiME HBL and EMR-Lipid.

The sample pH intensely influenced the efficiency with which the parabens and phenolic compounds were retained by the SPE sorbent. Tests were conducted over the pH range of 1–10 by passing 50 mL of aqueous solutions containing 500 ng/L of each target compound adjusted with diluted HCl or NaOH. It was found that the extraction yields of the analytes were the highest in the range between pH 3 and 5, and then diminished at higher values, so pH 4 was selected for further testing.

3.2. UHPLC-MS/MS Analysis

All of the analytes were infused and characterized by optimizing the multiple reaction monitoring (MRM) transitions for each one. As can be seen in Table 2, the precursor ion of each analyte was detected and two product transitions were optimized. The primary transition involved the most abundant product ion, which was used as the quantifying ion (Q), and a secondary ion that was used as the confirming ion (q). The procedure was applied with the software XCalibur 3.0.63 from Thermo Fisher Scientific (San José, CA, USA) and the data were processed with TraceFinder 3.2.5.12.0, also from Thermo Fisher Scientific. Table 2 shows the values of the MRM transitions, optimized collision energy, and lens voltage. The final detection conditions are described in Section 2.2.

3.3. Validation of the Method

The proposed methodology was validated by evaluating its performance with the most representative samples (milk, custard, yogurt, milkshakes, cheese, cream, butter, and margarine). The performance was assessed in terms of selectivity, linearity, matrix effects, precision, and recovery. These quality-related properties were determined by running calibration curves with a blank milk sample that was spiked with the analytes at variable concentrations and analysed as described in Section 2.4. Table 2 shows the most salient results. The method provided a linear response over the concentration range of 3.3–20,000 ng/kg and the correlation coefficients (r^2) exceeded 0.994 for all of the analytes. Limits of detection (LODs) and quantification (LOQs) were calculated in terms of the signal-to-noise ratio (S/N = 3 for LOD and S/N = 10 for LOQ). As can be seen from Table 2, LODs ranged from 1.0 to 20 ng/kg and LOQs coincided with the lower limit of each linear range. The limits of quantification of the proposed method varied between 3.3 and 65 ng/kg. These values were lower than the EU regulatory limits of the nethods included in Table 1 were lower than the EU regulatory limits [7,25,27].

The precision of the method was assessed in terms of the relative standard deviation (RSD), both intraday and interday, by spiking representative uncontaminated samples (n = 12) with the analytes at three different concentration levels (100, 500, and 2000 ng/kg). As can be seen from Table 3 and Table S1 (Supplementary Material) intraday precision was 2.03–8.91% and interday precision was 5.21–11.44%, so both were good and on a par with those for alternative methods [34,38].

	Compoundo								RSD) (%) ^a							
	Compounds	Μ	Milk		gurt	Cus	stard	Milk	shake	Cre	eam	Ch	eese	Bu	tter	Marg	arine
		Intraday	Interday	Intraday	Interday	Intraday	Interday	Intraday	Interday	Intraday	Interday	Intraday	Interday	Intraday	Interday	Intraday	Interday
	Methylparaben	5.11	8.22	5.31	10.21	5.05	10.16	5.06	9.03	2.55	9.75	4.26	9.26	4.01	7.12	2.36	7.42
s	Ethylparaben	4.26	9.33	4.43	11.15	5.06	8.03	3.83	7.23	2.13	7.98	5.62	7.32	5.25	8.92	4.53	8.54
en	Isopropylparaben	5.13	11.10	5.55	9.08	6.88	9.39	5.19	9.77	8.26	9.28	3.13	5.52	4.37	11.44	6.06	10.23
ab	Butylparaben	5.25	11.28	4.23	9.17	5.20	9.40	3.43	6.29	6.28	9.58	4.22	9.13	4.28	10.58	5.27	10.44
ar	Propylparaben	5.06	9.35	6.12	10.35	4.18	9.84	3.81	6.27	3.29	7.33	2.03	8.57	6.64	6.35	6.69	9.65
	Isobutylparaben	4.27	8.63	5.46	10.42	3.36	8.42	5.73	8.88	8.48	10.1	3.65	6.56	3.46	7.72	8.50	10.16
	Benzylparaben	3.19	9.46	5.88	9.44	4.32	8.03	6.52	9.06	5.55	9.08	4.42	9.04	4.56	6.96	4.04	6.34
	Pentylphenol	3.00	7.78	7.20	9.46	3.54	10.04	5.06	10.03	4.23	9.17	5.82	11.02	5.83	7.38	6.68	8.43
	4-Chloro-3-methylphenol	5.34	9.50	4.25	8.58	4.07	10.15	2.82	9.52	6.12	10.35	5.25	7.24	5.62	8.54	4.52	8.06
	4-tert-Octylphenol	4.16	8.44	2.13	7.40	6.48	9.35	8.07	10.44	5.46	10.42	6.52	8.33	5.75	8.94	3.06	8.03
	Phenol	3.42	7.65	5.52	9.06	4.27	10.13	5.08	8.32	7.26	11.03	3.21	7.75	3.06	6.25	2.82	9.82
	4-tert-Butylphenol	5.15	11.16	5.06	10.03	3.33	9.23	2.24	9.66	4.28	9.44	4.27	8.24	7.52	9.18	7.07	9.44
	Nonylphenol	7.06	10.23	4.82	9.82	4.66	7.55	3.55	9.63	5.35	9.38	4.17	9.12	5.54	10.07	4.01	7.10
	2-tert-Butylphenol	5.27	10.44	7.07	9.44	7.91	10.25	4.05	6.01	4.20	8.17	5.07	7.02	7.36	9.25	4.25	8.99
	Bisphenol F	4.69	9.65	4.08	8.22	4.12	8.02	3.06	8.03	2.93	7.89	4.56	6.45	3.12	6.32	5.37	11.04
	3,4-Dimehylphenol	7.50	10.16	4.24	9.66	6.32	9.33	4.75	9.95	3.05	8.53	3.04	6.34	4.23	6.64	6.28	10.38
ols	2,5-Dimethylphenol	3.56	7.76	3.55	9.13	7.84	10.07	9.41	7.91	5.27	10.26	4.68	9.43	4.45	6.87	7.64	6.33
en	4-Phenylphenol	5.23	10.31	6.13	9.21	6.86	8.25	8.02	4.12	6.38	10.39	4.75	5.24	3.32	7.02	4.46	8.70
ЪЧ	2-Phenylphenol	4.45	9.22	7.26	11.03	4.12	7.32	9.06	6.32	5.20	11.40	3.37	7.25	3.55	9.74	5.56	7.95
	Bisphenol S	5.27	10.26	4.28	9.44	4.13	9.64	3.12	7.32	4.01	7.11	2.61	5.24	4.32	6.64	4.18	9.84
	4-Chlorophenol	6.38	10.39	5.29	10.77	6.32	9.33	7.13	9.64	3.25	7.99	4.32	8.03	2.58	7.54	3.36	8.42
	4-Hexylphenol	5.20	11.40	3.48	5.29	4.37	9.44	6.06	10.23	5.38	11.54	5.54	10.04	4.78	10.96	5.34	9.50
	Bisphenol A	4.18	8.84	3.86	6.27	3.28	10.58	9.82	4.66	5.98	9.38	5.07	10.15	7.63	9.68	7.26	11.03
	Bisphenol B	2.36	7.42	4.73	8.86	4.08	8.22	4.12	8.02	4.24	9.66	3.48	9.35	4.96	8.02	4.28	9.44
	4-Heptylphenol	4.53	8.54	4.32	7.22	3.36	6.76	3.05	9.03	3.55	9.13	5.27	10.13	5.88	9.03	2.36	7.42
	Bisphenol Z	4.66	9.65	6.23	8.41	4.18	8.84	3.86	6.27	4.12	7.32	3.33	9.23	4.24	8.85	5.20	11.40
	Pentachlorophenol	4.48	9.27	3.44	5.58	2.61	4.24	4.32	6.64	3.18	8.84	8.91	11.30	5.85	7.23	6.38	10.39
	Triclosan	5.00	10.16	5.36	9.63	4.13	9.64	4.05	9.05	4.05	9.51	3.12	8.22	5.16	10.41	5.20	11.40

^a RSD: relative standard deviation (n = 12) for 100 ng/kg (p < 0.05).

Matrix effects (*ME*) can arise during chamber ionization (HESI) as a result of coextractives interfering with the electrospray ionization of the analytes and with an adverse impact on precision [32]. *ME* is a measure of the deviation of the slope of the matrix calibration curve from that of the standard curve [39]:

$$ME = \left[\left(\frac{slope \ of \ matrix - matched \ curve}{slope \ of \ in - solvent \ curve} \right) - 1 \right] \cdot 100$$

The effects were classified as negligible $[(0\%)-(\pm 10\%)]$, mild $[(\pm 10\%)-(\pm 20\%)]$, medium $[(\pm 20\%)-(\pm 50\%)]$, or strong (>±50%). As shown in Table S2 (Supplementary Material) and Figure 3, ME ranged from 1 to 10%, so the matrix effects were negligible. These results are better than the previously reported values (1–20%), which can be ascribed to our using the EMR-Lipid sorbent [22]. Therefore, using EMR-Lipid in combination with Oasis PRiME HLB efficiently reduced the matrix effect, as can be seen in Figure S2 (Supplementary Material).



Figure 3. Matrix effect of parabens and phenolic compounds in milk (**A**), yoghurt (**B**), cheese (**C**), and butter (**D**) encountered with EMR-lipid and Oasis PRiME HLB in combination. The blue squares in the picture correspond to the *MFs* of the different analytes.

The accuracy of the method was evaluated in terms of recovery, using samples spiked with a mixture of target compounds at three different concentration levels (100, 500, and 2000 ng/kg) that were processes as described in Section 2.4. The analyte recoveries thus obtained from the triplicate determinations (n = 3) are shown in Tables 4 and S1. As can be seen, they ranged from 91 to 105%. These recoveries were also better than those obtained in previous work without the EMR-Lipid sorbent (85–108%) [22]. Dualde et al. [25], Tuzimski and Szubartowski [27], and De Almeida Soares et al. [42] used alternative extraction and determination techniques for bisphenols and parabens in dairy products and obtained recoveries of 83–115%, 15–103%, and 74–112%, respectively.

	Compounds				Recove	ries (% \pm SD, 1	$t = 3)^{a}$		
	compounds	Milk	Yogurt	Cheese	Milkshake	Cream	Custard	Butter	Margarine
	Methylparaben	102 ± 8	100 ± 9	91 ± 8	95 ± 6	95 ± 7	95 ± 7	94 ± 6	91 ± 9
	Ethylparaben	103 ± 8	91 ± 6	95 ± 7	102 ± 7	102 ± 6	96 ± 9	98 ± 5	91 ± 9
sus	Isopropylparaben	95 ± 7	97 ± 5	95 ± 10	101 ± 6	101 ± 6	101 ± 9	104 ± 8	98 ± 10
rabe	Butylparaben	100 ± 9	99 ± 4	99 ± 10	101 ± 10	104 ± 7	93 ± 6	100 ± 5	99 ± 5
Pa	Propylparaben	104 ± 5	101 ± 6	95 ± 9	99 ± 9	99 ± 9	95 ± 6	104 ± 7	95 ± 9
	Isobutylparaben	96 ± 7	104 ± 7	105 ± 7	103 ± 9	99 ± 8	96 ± 11	99 ± 9	105 ± 7
	Benzylparaben	105 ± 6	105 ± 8	102 ± 8	91 ± 9	103 ± 6	95 ± 6	105 ± 5	100 ± 8
	Pentylphenol	103 ± 9	103 ± 5	101 ± 6	99 ± 9	98 ± 8	101 ± 6	100 ± 8	96 ± 11
	4-Chloro-3-methylphenol	96 ± 6	93 ± 6	100 ± 5	101 ± 10	100 ± 8	94 ± 7	103 ± 8	96 ± 11
	4-tert-Octylphenol	102 ± 7	92 ± 7	100 ± 5	96 ± 11	104 ± 8	95 ± 6	101 ± 8	106 ± 6
	Phenol	103 ± 9	97 ± 7	91 ± 8	100 ± 8	96 ± 6	95 ± 6	100 ± 7	103 ± 5
	4-tert-Butylphenol	101 ± 6	98 ± 8	97 ± 9	104 ± 8	94 ± 7	91 ± 9	103 ± 6	91 ± 8
	Nonylphenol	100 ± 8	91 ± 9	94 ± 7	102 ± 7	95 ± 6	96 ± 11	102 ± 7	104 ± 8
	2-tert-Butylphenol	104 ± 8	91 ± 9	95 ± 6	96 ± 6	100 ± 7	96 ± 6	103 ± 6	102 ± 7
	Bisphenol F	96 ± 11	99 ± 9	101 ± 6	96 ± 7	104 ± 7	96 ± 9	98 ± 8	97 ± 9
	3,4-Dimethylphenol	105 ± 7	99 ± 8	101 ± 10	101 ± 9	91 ± 8	98 ± 5	93 ± 7	94 ± 7
ols	2,5-Dimethylphenol	99 ± 8	103 ± 8	96 ± 9	98 ± 5	99 ± 5	104 ± 8	99 ± 9	98 ± 5
Pher	4-Phenylphenol	102 ± 7	102 ± 7	101 ± 8	103 ± 9	102 ± 6	98 ± 5	96 ± 6	99 ± 9
-	2-Phenylphenol	100 ± 7	96 ± 6	95 ± 6	103 ± 8	98 ± 5	99 ± 9	98 ± 8	96 ± 11
	Bisphenol S	98 ± 10	101 ± 9	106 ± 6	102 ± 7	95 ± 6	101 ± 8	94 ± 6	99 ± 9
	4-Chlorophenol	99 ± 5	100 ± 7	103 ± 5	98 ± 5	99 ± 9	95 ± 6	102 ± 9	94 ± 7
	4-Hexylphenol	102 ± 6	105 ± 8	96 ± 8	104 ± 8	96 ± 11	96 ± 6	102 ± 9	95 ± 6
	Bisphenol A	96 ± 9	101 ± 9	97 ± 7	98 ± 5	105 ± 7	98 ± 10	95 ± 7	100 ± 5
	Bisphenol B	101 ± 9	95 ± 7	97 ± 7	104 ± 8	98 ± 10	98 ± 5	105 ± 7	91 ± 8
	4-Heptylphenol	96 ± 9	102 ± 6	97 ± 6	97 ± 5	95 ± 10	104 ± 8	98 ± 5	103 ± 9
	Bisphenol Z	95 ± 6	95 ± 8	103 ± 6	99 ± 4	99 ± 10	97 ± 9	105 ± 10	96 ± 11
	Pentachlorophenol	99 ± 9	95 ± 8	99 ± 8	99 ± 8	103 ± 5	94 ± 7	97 ± 8	105 ± 7
	Triclosan	102 ± 5	97 ± 9	90 ± 8	102 ± 7	96 ± 8	104 ± 8	91 ± 9	104 ± 8

Table 4. Results obtained in the recovery studies of the different types of dairy product samples.

^a Percent recoveries ($\% \pm$ SD, *n* = 3) of analytes spiked to milk and dairy product samples for 100 ng/kg.

3.4. Application to Real Samples

Once the proposed the method was checked whether it performed well regarding sensitivity, recovery, precision, and matrix effects, it was applied to real samples. Thus, it was used to determine seven parabens and twenty-one phenolic compounds in fifteen samples of milk (skimmed, semi-skimmed and whole milk from cow; fresh cow milk; semi-skimmed goat milk; and semi-skimmed sheep milk) and sixteen samples of dairy products (milkshakes, yogurts, cheeses, butters, margarines, and creams).

All samples were analysed in triplicate and a blank was introduced after each three to discard potential contamination from other sources. Samples were chosen in terms of representativeness on the market and spanned various packaging materials including high-density polyethylene (HDPE), multilayer-cardboard (MC), polyethylene terephthalate (PET), polystyrene (PS), and glass.

As can be seen from Tables 5 and 6, the most frequently detected contaminants were BPA (68% of samples) and BPZ (39%). BPA was found at concentrations over the range 33–580 ng/kg, especially in whole cow milk and milkshakes held in MC, PET, or HDPE. In contrast, it was completely absent from yogurt samples held in glass. De Almeida Soares et al. previously detected BPA at concentrations of 77,600–150,800 ng/L in milk [38], and Lv et al. found it at levels from 800 to 128,700 g/kg in dairy products [37]. Azzouz et al. detected BPA at concentrations from 1400 to 2900 ng/L in human milk, thus demonstrating that contamination from food and the environment can reach human bodies [30].

	Skimmed Cow's Milk ^b						Semi-Skimmed Cow's Milk			Semi-Skimmed Goat'S Milk		Semi-Skimmed Sheep's Milk	Whole Cow's Milk			
	Sample (Packaging Material)	1 (MC) °	2 (MC)	3 (MC)	4 (HDPE)	1 (MC)	2 (MC)	3 (HDPE)	4 (PET)	1 (MC)	2 (MC)	1 (MC)	1 (MC)	2 ° (MC)	3 (PET)	4 (HDPE)
	Nonylphenol	_ d	-	-	-	-	-	130 ± 10	34 ± 3	210 ± 20	-	53 ± 4	-	-	-	-
Ś	4-tert-butylphenol	170 ± 10	-	-	-	81 ± 6	220 ± 20	110 ± 10	21 ± 1	220 ± 20	-	-	-	-	-	-
ed	Ethylparaben	80 ± 6	63 ± 3	47 ± 4	-	33 ± 2	-	-	-	470 ± 40	300 ± 30	210 ± 20	300 ± 20	-	-	-
on	4-Chlorophenol	160 ± 10	110 ± 10	60 ± 5	-	50 ± 4	-	-	-	130 ± 10	-	-	-	-	-	-
np et	Propylparaben	-	-	-	-	-	10 ± 1	-	9 ± 0.8	33 ± 3	24 ± 2	23 ± 2	30 ± 2	-	-	-
P 5	Benzyl-paraben	-	-	-	-	31 ± 3	24 ± 2	-	33 ± 1	98 ± 9	90 ± 6	86 ± 6	130 ± 10	-	-	-
•	Bisphenol A	130 ± 10	33 ± 3	-	-	67 ± 5	-	-	-	64 ± 5	99 ± 8	170 ± 10	580 ± 50	260 ± 20	190 ± 10	290 ± 20
	Bisphenol Z	33 ± 3	44 ± 4	41 ± 4	25 ± 2	35 ± 3	34 ± 3	24 ± 2	36 ± 3	57 ± 5	47 ± 4	45 ± 4	42 ± 3	-	-	-
	Pentachlorophenol	-	-	-	-	-	40 ± 3	35 ± 3	55 ± 5	40 ± 2	-	76 ± 6	-	-	-	-

Table 5. Parabens and phenolic compound concentrations (mean \pm standard deviation, ng/kg, n = 3) detected in different milk^a.

^a The compounds not listed were not detected. ^b (fat content,%)/(protein content,%): skimmed cow's milk: (0.05–0.3%)/(3.2–3.4%); semi-skimmed cow's milk (0.05–0.3%)/(3.2–3.4%); semi-skimmed sheep 's milk (1.6%)/(5.4%); whole cow's milk (3.6%)/(3.0–3.2%); ^c HDPE: High density polyethylene; MC: multilayer-cardboard; PET: polyethylene terephthalate; PS: polystyrene. ^d -: not detected. ^e Fresh whole cow's milk.

Table 6. Parabens and phenolic compound concentrations (mean \pm standard deviation, ng/kg, n = 3) detected in dairy products ^a.

		Dairy Products															
				Yogurt							D (/		Managering		C		
			MIIKSNAKe		Sheep	Cow	Cow	Sheep	-	Cheese		Butter		wargarine		Cream	
ted unds	Sample (packaging material)	1 (PET)	2 (HDPE)	3 (MC)	1 (Glass)	2 (Glass)	3 (PS)	4 (Glass)	1 (PS)	2 (PS)	3 (PS)	1 (PS)	2 (PS)	1 (PS)	2 (PS)	1 (MC)	2 (MC)
etec	Nonylphenol	_ d	-	-	43 ± 3	-	-	-	-	-	-	-	98 ± 8	99 ± 7	49 ± 4	-	-
D 10	Bisphenol A	360 ± 20	420 ± 40	370 ± 30	-	-	58 ± 5	-	-	150 ± 10	-	140 ± 10	50 ± 5	170 ± 10	160 ± 10	260 ± 20	450 ± 40

^a The compounds not listed were not detected. ^b (fat content,%)/(protein content,%): milkshake (0.8-0.9%)/(1.6-3%); yogurt (1.9-6.4%)/(3.1-5.5%); cheese (14-15%)/(10-16%); butter (80-82%)/(0.4-0.6%); margarine (60-70%)/(0.3%); cream (18-38%) (1.8–2.5\%). ^c HDPE: High density polyethylene; MC: multilayer-cardboard; PET: polyethylene terephthalate; PS: polystyrene. ^d -: not detected.

BPZ was the second most detected contaminant at 24 to 57 ng/kg. However, as displayed in Tables 5 and 6, it was only detected in milk samples, with its peak concentrations in the samples held in multilayer cardboard (MC) and its lowest levels in those held in HDPE. No BPZ analogues, such as BPF, BPS, or BPB were detected. In contrast, De Almeida Soares et al. [42] and Di Marco Pisciottato et al. [38] found BPS at concentrations of 60,000 ng/L and BPF at levels of 500–8700 ng/L, respectively, in milk.

Nonylphenol is one other frequently monitored contaminant. In 2022, it was included into the candidate contaminants for the regulatory determination list [20]. This phenolic compound was found in 26% of our samples, at concentrations over the range of 34–210 ng/kg in milk packed in HDPE, PET, and MC, and in both butter and margarine packed in PS. As displayed in Table 1, Lv et al. detected nonylphenol at concentrations of 3500–36,700 ng/kg in dairy products [37] and Shao et al. found it at 100 ng/kg in milk samples [26].

4-*tert*-Butylphenol was spotted in 19% of the milk samples (concentrations of 21–220 ng/kg), and 4-chlorophenol (50–160 ng/kg) and pentachlorophenol (35–76 ng/kg) both in 16%. As shown in Table 5, these contaminants were present in the majority of samples of milk packed in MC. Zhang et al. previously detected 4-*tert*-butylphenol at 850 ng/kg in milk [23].

Parabens are employed as food preservatives for their anti-microbial and anti-fungal properties. As can be seen in Table 3, only three of the parabens studied were detected in the samples. The following were only found in milk, regardless of the packaging material: ethylparaben (33–470 ng/kg), propylparaben (9–33 ng/kg), and benzylparaben (24–130 ng/kg). Dualde et al. [25] and Azzouz et al. [30] detected parabens at concentrations of 130–4050 ng/L and 15–8100 ng/L, respectively, in breast milk.

4. Conclusions

A new method for the simultaneous determination of 21 phenol compounds and 7 parabens in various types of dairy products was developed. The method combines conventional LLE extraction with SPE in a closed semi-automated system furnished with two different SPE sorbents (ERM-lipid and Oasis PRiME HLB) to avoid unwanted matrix effects, and increase analyte recoveries and expedite sample treatment. This substantially reduced the consumption of samples and solvents, as well as the analysis time, relative to the existing methods for the same purpose. Thus, the proposed method speeds up the sample treatment by reducing the time of this step to 15 min compared with other methods that need more than 40 min [4,24,25]. The optimized UHPLC-MS/MS method allowed the target analytes to be detected and quantified with good precision (RSD 2-11.4%), low matrix effects (0–10%), and high recoveries (91–105%). These results are quite good if one considers the high complexity of the sample matrices. A total of 31 samples (15 of milk and 16 of dairy products) were analysed. The most detected phenolic compounds were BPA (33-580 ng/kg) and BPZ (24-57 ng/kg). In the case of BPA, these concentration levels did not exceed the legal limit set by the European Union. Nonylphenol was found in 26% of the samples, at concentrations up to 210 ng/kg. The most detected paraben was ethylparaben, which was found at levels up to 470 ng/kg. Glass packaging was found to hold the least contaminated samples and plastics the most contaminated ones (particularly milk).

Supplementary Materials: The following supporting information can be downloaded at: https:// www.mdpi.com/article/10.3390/foods12152909/s1, Figure S1. Continuous-flow system for the solidphase with three injection valves and two serially columns; Figure S2. Extract ion chromatograms for four analytes (benzylparaben, bisphenol A, triclosan and 4-heptylphenol). In one chromatogram, the analyte is spiked in the solvent, and in the other the analyte is spiked in a milk sample; Table S1. Results obtained for precision and recovery studies for the different types of dairy product samples; Table S2. Results obtained in the matrix effect studies of the different types of dairy product samples. Author Contributions: Conceptualization, L.P.C.; methodology, L.P.C., A.J.R. and E.B.; investigation, L.P.C., A.J.R. and E.B.; Validation, L.P.C., A.J.R. and E.B.; writhing—original draft, L.P.C.; writing—review and editing, A.J.R. and E.B.; funding acquisition, E.B.; supervision, E.B. All authors have read and agreed to the published version of the manuscript.

Funding: The authors received funding from the Consejería de Economía, Conocimiento, Empresas y Universidad, Regional Government of Andalucía, Spain (Project Ref. PY2018-1211), and partially supported by EU FEDER funds, and the Research Program of the University of Jaen (Plan 2019–2020).

Data Availability Statement: The data used to support the findings of this study can be made available by the corresponding author upon request.

Acknowledgments: Technical and human support provided by CICT of University of Jaén (UJA, MINECO, Junta de Andalucía, FEDER) is gratefully acknowledged.

Conflicts of Interest: The authors declare no conflict of interest.

References

- Ministry of Agriculture, Fisheries and Food. Milk Month by Month. Catalog of Publications of the General State. 2021. Available online: https://www.mapa.gob.es/es/alimentacion/temas/consumo-tendencias/la-leche-mes-a-mes-agosto-2021_tcm30-58 1083.pdf (accessed on 30 June 2023).
- 2. Vale, F.; Sousa, C.A.; Sousa, H.; Santos, L.; Simoes, M. Parabens as emerging contaminants: Environmental persistence, current practices and treatment processes. *J. Clean. Produ.* 2022, 347, 131244. [CrossRef]
- Pan, Y.; Liu, X.; Liu, J.; Wang, J.; Liu, J.; Gao, Y.; Ma, N. Determination of organophosphorus in dairy products by graphitic carbon nitride combined molecularly imprinted microspheres with ultra performance liquid chromatography. *Food Chem.* 2022, 15, 100424. [CrossRef] [PubMed]
- 4. Wang, J.; Zhao, W.; Guo, W.; Li, Y.; Jiang, R.; Li, H.; Wang, S.; Li, Z. Simultaneous screening and analysis of 155 veterinary drugs in livestock foods using ultra-high performance liquid chromatography tandem quadrupole linear-ion-trap mass spectrometry. *Food Chem.* **2022**, *393*, 133260. [CrossRef] [PubMed]
- Deng, Z.H.; Li, N.; Jiang, H.L.; Lin, J.M.; Zhao, R.S. Pretreatmetechniques and analytical methods for phenolic endocrine disrupting chemicals in food and environmental samples. *TrAC—Trends Anal. Chem.* 2019, 119, 115592. [CrossRef]
- 6. Gu, L.; Wang, Z.; Pan, Y.; Wang, H.; Sun, L.; Liu, L.; Zhong, W.J.; Zhang, W.B. Associations between mixed urinary phenols and parabens metabolites and bone mineral density: Four statistical models. *Chemosphere* **2023**, *31*, 137065. [CrossRef]
- Herrero, L.; Quintanilla-López, J.E.; Fernández, M.A.; Gómara, B. Plasticisers and preservatives in commercial milk products: A comprehensive study on packages used in the Spanish market. *Food Chem.* 2021, 338, 128031. [CrossRef]
- Iribarne-Durán, L.M.; Serrano, L.; Peinado, F.M.; Peña-Caballero, M.; Hurtado, J.A.; Vela-Soria, F.; Fernández, M.F.; Freire, C.; Artacho-Cordóng, F.; Olea, N. Biomonitoring bisphenols, parabens, and benzophenones in breast milk from a human milk bank in Southern Spain. *Sci. Total Environ.* 2022, 830, 154737. [CrossRef]
- 9. Farajzadeh, M.A.; Nemati, M.; Altunay, N.; Tuzen, M.; Kaya, S.; Kheradmand, F. Experimental and density functional theory studies during a new solid phase extraction of phenolic compounds from wastewater samples prior to GC–MS determination. *Microchem. J.* 2022, *177*, 107291. [CrossRef]
- Marazuela, M.D.; Klaiber, M.; Moreno-Gordaliza, E.; Barata, A.; Gómez-Gómez, M.M. Safety assessment of commercial antimicrobial food packaging: Triclosan and microplastics, a closer look. *Food Packag. Shelf Life* 2022, *31*, 100780. [CrossRef]
- Caballero-Casero, N.; Lunar, L.; Rubio, S. Analytical methods for the determination of mixtures of bisphenols and derivatives in human and environmental exposure sources and biological fluids. A review. *Anal. Chim. Acta* 2016, 908, 22–53. [CrossRef]
- Zhou, Q.; Lei, M.; Li, J.; Zhao, K.; Liu, Y. Sensitive determination of bisphenol A, 4-nonylphenol and 4-octylphenol by magnetic solid phase extraction with Fe@MgAl-LDH magnetic nanoparticles from environmental water samples. *Sep. Purif. Technol.* 2017, 182, 78–86. [CrossRef]
- Adjei, J.K.; Ofori, A.; Megbenu, H.K.; Ahenguah, T.; Boateng, A.K.; Adjei, G.; Bentum, J.K.; Essumang, D.K. Health risk and source assessment of semi-volatile phenols, p-chloroaniline and plasticizers in plastic packaged (sachet) drinking water. *Sci. Total Environ.* 2021, 797, 149008. [CrossRef]
- EFSA. Bisphenol A: EFSA Draft Opinion Proposes Lowering the Tolerable Daily Intake. European Food Safety Agency. 2021. Available online: https://www.efsa.europa.eu/en/news/bisphenol-efsa-draft-opinion-proposes-lowering-tolerable-daily-intake (accessed on 30 June 2023).
- 15. European Commission. Commission Regulation (EU) 2018/213 of 12 February 2018 on the use of bisphenol A in varnishes and coatings intended to come into contact with food and amending Regulation (EU) No 10/2011 as regards the use of that substance in plastic food contact. *Off. J. Eur. Community* **2018**, *L*41, 6–12.
- 16. Tuzimski, T.; Szubartowski, S. Application of d-SPE before spe and HPLC-FLD to analyze bisphenols in human breast milk samples. *Molecules* **2021**, *26*, 4930. [CrossRef]
- 17. Santonicola, S.; Albrizio, S.; Ferrante, M.C.; Mercogliano, R. Study on bisphenol F, a bisphenol A analogue, at a dairy company: Health hazard and risk assessment. *Food Chem. Toxicol.* **2021**, *154*, 112334. [CrossRef] [PubMed]

- European Commission. Regulation (EU) No 10/2011 of 14 January 2011 on plastic materials and articles intended to come into contact with food. Off. J. Eur. Community 2011, L12, 1–89.
- USEPA. Appendix A to 40 CFR, Part 423—126 Priority Pollutants; EPA: Washington, DC, USA, 2022. Available online: https://www.ecfr.gov/current/title-40/chapter-I/subchapter-N/part-423/appendix-Appendix%20A%20to%20Part%20423 (accessed on 30 June 2023).
- USEPA. Contaminant Candidate List (CCL) and Regulatory Determination. Chemical Contaminants—CCL 4. Final CCL 4 Chemical Contaminants. 2022. Available online: https://www.epa.gov/ccl/chemical-contaminants-ccl-4 (accessed on 30 June 2023).
- González, N.; Cunha, S.C.; Ferreira, R.; Fernandes, J.O.; Marqués, M.; Nadal, M.; Domingo, J.L. Concentrations of nine bisphenol analogues in food purchased from Catalonia (Spain): Comparison of canned and non-canned foodstuffs. *Food Chem. Toxicol.* 2020, 136, 110992. [CrossRef]
- 22. Palacios Colón, L.; Rascón, A.J.; Ballesteros, E. Simultaneous determination of phenolic pollutants in dairy products held in various types of packaging by gas chromatography—mass spectrometry. *Food Control* **2023**, *146*, 109564. [CrossRef]
- Zhang, J.; Dang, X.; Dai, J.; Hu, Y.; Chen, H. Simultaneous detection of eight phenols in food contact materials after electrochemical assistance solid-phase microextraction based on amino functionalized carbon nanotube/polypyrrole composite. *Anal. Chim. Acta* 2021, 1183, 338981. [CrossRef]
- Xiong, L.; Yan, P.; Chu, M.; Gao, Y.Q.; Li, W.H.; Yang, X.L. A rapid and simple HPLC–FLD screening method with QuEChERS as the sample treatment for the simultaneous monitoring of nine bisphenols in milk. *Food Chem.* 2018, 244, 371–377. [CrossRef] [PubMed]
- Dualde, P.; Pardo, O.F.; Fernández, S.; Pastor, A.; Yusà, V. Determination of four parabens and bisphenols A, F and S in human breast milk using QuEChERS and liquid chromatography coupled to mass spectrometry. J. Chromatogr. B 2019, 1114–1115, 154–166. [CrossRef]
- 26. Shao, B.; Han, H.; Tu, X.; Huang, L. Analysis of alkylphenol and bisphenol A in eggs and milk by matrix solid phase dispersion extraction and liquid chromatography with tandem mass spectrometry. J. Chromatogr. B 2007, 850, 412–416. [CrossRef] [PubMed]
- Tuzimski, T.; Szubartowski, S. Method development for selected bisphenols analysis in sweetened condensed milk from a can and breast milk samples by HPLC–DAD and HPLC-QQQ-MS: Comparison of sorbents (Z-Sep, Z-Sep plus, PSA, C18, Chitin and EMR-lipid) for clean-up of QuEChERS extract. *Molecules* 2019, 24, 2093. [CrossRef]
- Ye, X.; Bishop, A.M.; Needham, L.L.; Calafat, A.M. Automated on-line column-switching HPLC-MS/MS method with peak focusing for measuring parabens, triclosan, and other environmental phenols in human milk. *Anal. Chim. Acta* 2008, 622, 150–156. [CrossRef] [PubMed]
- Gao, Y.; Wang, S.; Zhang, N.; Xu, X.; Bao, T. Novel solid-phase extraction filter based on a zirconium meta-organic framework for determination of non-steroidal anti-inflammatory drugs residues. J. Chromatogr. A 2021, 1652, 462349. [CrossRef]
- 30. Azzouz, A.; Rascón, A.J.; Ballesteros, E. Simultaneous determination of parabens, alkylphenols, phenylphenols, bisphenol A and triclosan in human urine, blood and breast milk by continuous solid-phase extraction and gas chromatography-mass spectrometry. *J. Pharm. Biomed. Anal.* **2016**, *119*, 16–26. [CrossRef]
- Castilla-Fernández, D.; Moreno-González, D.; Beneito-Cambra, M.; Molina-Díaz, A. Critical assessment of two sample treatment methods for multiresidue determination of veterinary drugs in milk by UHPLC-MS/MS. *Anal. Bioanal. Chem.* 2019, 411, 1433–1442. [CrossRef]
- Yin, Z.; Yuan, B.; Lyu, W.; Huang, Q.; Simon, J.E.; Wu, Q. Method development and validation for analysis of phenolic compounds in fatty complex matrices using enhanced matrix removal (EMR) lipid cleanup and UHPLC-QqQ-MS/MS. *Food Chem.* 2022, 373, 131096. [CrossRef]
- Zhao, L.; Lukas, D. Multiresidue Analysis of Veterinary Drugs in Bovine Liver by LC/MS/MS Agilent Bond Elut Enhanced Matrix Removal—Lipid, Agilent Technologies Application Note. 2015. Available online: https://www.agilent.com/cs/library/ applications/5991-6096EN.pdf (accessed on 30 June 2023).
- Reyes-Gallardo, E.M.; Lucena, R.; Cárdenas, S.; Valcárcel, M. Dispersive micro-solid phase extraction of bisphenol A from milk using magnetic nylon 6 composite and its final determination by HPLC-UV. *Microchem. J.* 2016, 124, 751–756. [CrossRef]
- 35. Qiao, L.; Sun, R.; Tao, Y.; Yan, Y. New low viscous hydrophobic deep eutectic solvents for the ultrasound-assisted dispersive liquid-liquid microextraction of endocrine-disrupting phenols in water, milk and beverage. *J. Chromatogr. A* 2022, 1662, 462728. [CrossRef]
- Matuszewski, B.K.; Constanzer, M.L.; Chavez-Eng, C.M. Strategies for the assessment of matrix effect in quantitative bioanalytical methods based on HPLC-MS/MS. *Anal. Chem.* 2003, 75, 3019–3030. [CrossRef] [PubMed]
- Lv, T.; Zhao, X.E.; Zhu, S.; Qu, F.; Song, C.; You, J.; Suo, Y. Determination of bisphenol A, 4-octylphenol, and 4-nonylphenol in soft drinks and dairy products by ultrasound-assisted dispersive liquid-liquid microextraction combined with derivatization and high-performance liquid chromatography with fluorescence detection. J. Sep. Sci. 2014, 37, 2757–2763. [CrossRef]
- 38. Di Marco Pisciottano, I.; Guadagnuolo, G.; Busico, F.; Alessandroni, L.; Neri, B.; Vecchio, D.; Di Vuolo, G.; Cappelli, G.; Martucciello, A.; Gallo, P. Determination of 20 endocrine-disrupting compounds in the buffalo milk production chain and commercial bovine milk by UHPLC–MS/MS and HPLC–FLD. *Animals* **2022**, *12*, 410. [CrossRef] [PubMed]

- Czarczynska-Goslinska, B.; Grzeskowiak, T.; Frankowski, R.; Lulek, J.; Pieczak, J.; Zgoła-Grzeskowiak, A. Determination of bisphenols and parabens in breast milk and dietary risk assessment for Polish breastfed infants. *J. Food Compost. Anal.* 2021, *98*, 103839. [CrossRef]
- 40. Muñoz Hornillos, M.; Yoldi Bienzobas, G. Leche y derivados. In *Alimentos Composición y Propiedades*; Astiasarán, I., Anchía, J.A., Hernández, M., Eds.; McGraw-Hill–Interamericana de España, S.A.U: Madrid, Spain, 2003; pp. 69–108.
- 41. Liao, C.; Kannan, K. Concentrations and profiles of bisphenol A and other bisphenol analogues in foodstuffs from the United States and their implications for human exposure. *J. Agric. Food Chem.* **2013**, *61*, 4655–4662. [CrossRef] [PubMed]
- 42. De Almeida Soares, D.; Pereira, I.; Sousa, J.C.P.; Bernardo, R.A.; Simas, R.C.; Vaz, B.G.; Chaves, A.R. Bisphenol determination in UHT milk and packaging by paper spray ionization mass spectrometry. *Food Chem.* **2023**, *400*, 134014. [CrossRef]

Disclaimer/Publisher's Note: The statements, opinions and data contained in all publications are solely those of the individual author(s) and contributor(s) and not of MDPI and/or the editor(s). MDPI and/or the editor(s) disclaim responsibility for any injury to people or property resulting from any ideas, methods, instructions or products referred to in the content.