

Supplementary Material

Determination of Parabens and Phenolic Compounds in Dairy Products through the Use of a Two-Step Continuous SPE System Including an Enhanced Matrix Removal Sorbent in Combination with UHPLC–MS/MS

Laura Palacios Colón, Andrés J. Rascón and Evaristo Ballesteros *

Department of Physical and Analytical Chemistry, E.P.S of Linares, University of Jaén,
23700 Linares, Jaén, Spain

Table S1. Results obtained for precision and recovery studies for the different types of dairy product samples.

		RSD (%) /Recovery (%) ^a															
		Milk		Yogurt		Custard		Milkshake		Cream		Cheese		Butter		Margarine	
		500 ng/kg	2000 ng/kg	500 ng/kg	2000 ng/kg	500 ng/kg	2000 ng/kg	500 ng/kg	2000 ng/kg	500 ng/kg	2000 ng/kg	500 ng/kg	2000 ng/kg	500 ng/kg	2000 ng/kg	500 ng/kg	2000 ng/kg
Parabens	Methylparaben	6.10 ^a	5.28	5.33	10.25	6.46	7.21	7.50	8.67	8.25	9.75	7.25	9.28	6.09	7.19	6.05	7.45
		99 ± 10 ^a	101 ± 10	104 ± 7	101 ± 9	95 ± 7	95 ± 7	98 ± 10	101 ± 9	105 ± 9	102 ± 7	96 ± 11	93 ± 7	95 ± 6	102 ± 9	101 ± 6	106 ± 6
	Ethylparaben	9.20	7.16	6.30	11.40	5.39	8.08	7.86	7.27	5.85	7.98	5.69	7.35	6.44	8.91	5.52	8.51
		101 ± 6	99 ± 10	91 ± 8	103 ± 8	105 ± 9	98 ± 5	103 ± 9	92 ± 6	91 ± 8	103 ± 9	101 ± 10	96 ± 8	101 ± 10	97 ± 6	101 ± 6	103 ± 7
	Isopropylparaben	5.28	7.35	7.55	9.87	6.80	9.39	5.33	9.73	8.47	9.50	9.05	5.52	8.34	11.46	6.00	10.37
		103 ± 9	103 ± 10	93 ± 6	97 ± 7	98 ± 10	104 ± 8	95 ± 7	95 ± 7	94 ± 6	91 ± 9	97 ± 6	103 ± 9	103 ± 9	96 ± 9	102 ± 8	95 ± 9
	Butylparaben	5.36	7.83	8.05	7.84	7.21	9.47	6.46	6.27	5.48	9.58	9.38	9.17	5.69	10.81	5.49	10.44
		98 ± 5	99 ± 9	104 ± 8	99 ± 10	98 ± 5	101 ± 6	99 ± 9	103 ± 9	95 ± 6	99 ± 9	104 ± 8	99 ± 9	102 ± 7	104 ± 8	99 ± 6	92 ± 7
	Propylparaben	8.01	9.56	6.21	10.57	7.18	9.87	8.80	6.20	7.87	7.38	11.07	8.54	6.64	6.58	6.62	9.50
		104 ± 8	98 ± 10	104 ± 8	95 ± 6	101 ± 8	106 ± 6	102 ± 9	95 ± 6	104 ± 9	103 ± 9	96 ± 9	98 ± 5	96 ± 6	95 ± 6	102 ± 9	102 ± 7
	Isobutylparaben	8.29	7.30	9.66	10.22	7.37	8.42	5.77	8.81	8.48	10.14	7.10	6.54	6.49	7.26	8.59	10.60
		102 ± 7	96 ± 8	94 ± 6	99 ± 9	104 ± 7	101 ± 6	97 ± 5	98 ± 9	99 ± 6	99 ± 10	104 ± 8	101 ± 6	101 ± 8	103 ± 8	99 ± 4	103 ± 10
	Benzylparaben	8.19	8.69	5.95	7.15	8.04	6.82	7.27	8.09	6.82	9.08	5.71	9.62	4.95	7.64	7.04	6.45
		99 ± 7	104 ± 7	93 ± 6	102 ± 6	97 ± 6	98 ± 7	104 ± 9	101 ± 6	101 ± 9	104 ± 8	98 ± 10	104 ± 8	99 ± 10	100 ± 6	99 ± 8	98 ± 7


Table S1. Continued

	RSD (%) / Recovery (%) ^a																
	Milk		Yogurt		Custard		Milkshake		Cream		Cheese		Butter		Margarine		
	500 ng/kg	2000 ng/kg	500 ng/kg	2000 ng/kg	500 ng/kg	2000 ng/kg	500 ng/kg	2000 ng/kg	500 ng/kg	2000 ng/kg	500 ng/kg	2000 ng/kg	500 ng/kg	2000 ng/kg	500 ng/kg	2000 ng/kg	
Phenols	Pentylphenol	6.49 ^a	7.95	7.20	8.77	6.55	10.40	5.69	5.78	7.29	8.16	8.57	8.52	7.83	7.88	6.66	8.30
		98 ± 6 ^a	99 ± 5	104 ± 8	97 ± 5	101 ± 6	99 ± 8	99 ± 8	97 ± 6	94 ± 7	104 ± 8	105 ± 10	98 ± 9	95 ± 7	100 ± 5	99 ± 8	105 ± 7
	4-Chloro-3-methylphenol	5.29	9.56	5.20	5.27	10.67	4.65	9.38	9.55	6.10	10.54	8.48	5.34	10.15	6.19	9.15	6.89
		104 ± 8	99 ± 9	105 ± 7	98 ± 9	98 ± 10	101 ± 9	91 ± 8	101 ± 10	93 ± 6	95 ± 7	92 ± 7	103 ± 5	99 ± 8	94 ± 7	103 ± 9	103 ± 9
	4- <i>tert</i> -Octylphenol	5.39	8.49	5.65	4.34	7.46	9.50	4.20	10.47	5.40	10.21	8.45	8.33	5.75	8.49	3.59	6.57
		99 ± 9	95 ± 8	99 ± 8	99 ± 8	103 ± 5	104 ± 8	102 ± 7	92 ± 6	94 ± 6	99 ± 9	94 ± 7	93 ± 6	103 ± 9	105 ± 9	98 ± 8	104 ± 9
	Phenol	8.08	8.37	5.95	7.91	9.57	10.34	6.14	8.21	7.24	11.30	6.52	7.53	3.08	6.56	8.88	3.25
		104 ± 9	101 ± 6	101 ± 9	104 ± 8	98 ± 10	104 ± 8	103 ± 9	99 ± 6	98 ± 10	92 ± 7	105 ± 8	105 ± 7	98 ± 10	95 ± 7	100 ± 5	103 ± 5
	4- <i>tert</i> -Butylphenol	5.59	11.1	5.44	6.36	7.39	9.34	5.29	9.65	6.27	9.40	9.86	8.40	7.51	9.82	7.01	9.40
		96 ± 6	103 ± 6	102 ± 7	100 ± 8	99 ± 8	98 ± 5	101 ± 6	99 ± 9	103 ± 9	96 ± 9	98 ± 8	97 ± 9	96 ± 9	104 ± 8	102 ± 9	104 ± 8
	Nonylphenol	7.64	10.31	5.84	9.24	6.65	7.54	5.55	9.37	5.32	9.81	10.77	9.25	5.58	10.74	6.08	7.17
		98 ± 10	101 ± 9	96 ± 11	101 ± 9	101 ± 6	99 ± 9	103 ± 9	103 ± 9	101 ± 10	104 ± 7	93 ± 6	102 ± 6	97 ± 6	101 ± 9	95 ± 6	91 ± 8
	2- <i>tert</i> -Butylphenol	5.70	10.44	7.02	9.41	7.57	10.51	7.08	6.17	6.22	9.61	6.38	9.34	5.78	9.50	7.24	8.97
		99 ± 8	102 ± 6	102 ± 7	99 ± 9	98 ± 10	95 ± 7	100 ± 5	96 ± 9	98 ± 5	99 ± 9	104 ± 8	95 ± 6	101 ± 8	106 ± 6	99 ± 9	101 ± 9
	Bisphenol F	5.99	9.55	6.04	8.27	9.41	8.24	5.62	8.3	5.89	7.91	6.50	6.57	6.65	6.24	5.34	11.31
		101 ± 6	94 ± 5	91 ± 8	101 ± 10	93 ± 6	93 ± 7	104 ± 8	103 ± 5	100 ± 5	103 ± 9	94 ± 7	103 ± 9	101 ± 10	96 ± 11	102 ± 7	104 ± 8
	3,4-Dimethylphenol	7.00	10.64	4.29	9.61	8.25	9.35	5.77	9.58	5.56	8.36	6.52	6.40	8.28	6.49	6.24	10.81
		100 ± 5	101 ± 9	95 ± 7	98 ± 10	97 ± 5	104 ± 8	97 ± 5	101 ± 6	93 ± 6	104 ± 8	98 ± 10	98 ± 5	99 ± 8	105 ± 7	103 ± 9	103 ± 9
	2,5-Dimethylphenol	6.69	7.71	5.51	9.32	9.09	10.71	9.81	7.16	5.44	10.64	5.67	9.38	6.43	6.78	7.60	6.38
		100 ± 8	99 ± 8	98 ± 5	95 ± 7	101 ± 6	106 ± 6	103 ± 9	100 ± 8	105 ± 7	91 ± 6	103 ± 6	103 ± 9	97 ± 8	104 ± 8	91 ± 9	95 ± 6
	4-Phenylphenol	5.33	10.13	6.15	9.19	11.30	8.56	8.03	7.23	6.32	10.98	6.76	5.45	5.29	7.25	5.47	8.09
		99 ± 5	93 ± 6	100 ± 5	94 ± 7	100 ± 8	91 ± 8	102 ± 7	96 ± 8	97 ± 6	97 ± 5	95 ± 10	100 ± 5	93 ± 6	100 ± 5	93 ± 6	103 ± 8
	2-Phenylphenol	6.32	9.25	7.25	11.30	6.17	7.26	9.03	6.24	5.25	11.01	5.68	7.55	8.98	9.41	5.55	7.56
		96 ± 6	93 ± 6	102 ± 7	104 ± 6	101 ± 9	104 ± 8	98 ± 5	103 ± 9	99 ± 4	101 ± 10	100 ± 5	93 ± 6	100 ± 5	91 ± 9	101 ± 6	104 ± 7
	Bisphenol S	4.35	10.61	8.26	9.41	6.19	9.45	7.58	6.19	8.94	7.98	8.32	5.44	6.34	6.47	7.18	9.41
		103 ± 9	104 ± 8	103 ± 9	105 ± 8	102 ± 8	99 ± 8	101 ± 10	101 ± 9	91 ± 9	99 ± 9	95 ± 8	99 ± 8	105 ± 9	95 ± 8	99 ± 8	99 ± 8
	4-Chlorophenol	6.28	10.90	5.24	10.75	6.35	9.30	7.14	6.30	5.99	6.35	7.57	9.31	5.56	7.41	6.55	8.20
		102 ± 7	98 ± 5	102 ± 6	104 ± 7	105 ± 7	96 ± 8	96 ± 6	102 ± 9	98 ± 5	105 ± 7	99 ± 9	104 ± 8	99 ± 9	98 ± 5	101 ± 10	101 ± 9
	4-Hexylphenol	9.09	11.09	6.46	5.28	8.33	9.41	6.08	6.39	5.31	9.35	6.15	10.41	7.74	10.61	5.35	9.08
		103 ± 9	102 ± 6	98 ± 5	99 ± 9	91 ± 9	94 ± 7	105 ± 7	98 ± 10	105 ± 7	99 ± 10	95 ± 7	100 ± 5	100 ± 7	93 ± 6	100 ± 5	94 ± 7
	Bisphenol A	5.34	8.42	6.86	6.28	7.56	10.81	9.83	9.83	5.91	9.89	7.07	10.50	7.64	9.60	7.21	11.38
		98 ± 8	99 ± 8	101 ± 10	101 ± 9	91 ± 9	98 ± 5	93 ± 7	94 ± 7	97 ± 7	98 ± 5	100 ± 5	93 ± 6	100 ± 5	96 ± 6	103 ± 6	102 ± 7
	Bisphenol B	6.63	7.20	5.68	5.55	10.45	8.27	6.35	8.27	5.42	9.60	6.70	9.51	5.94	8.20	6.25	9.40
		103 ± 5	101 ± 6	99 ± 9	98 ± 8	103 ± 7	103 ± 8	96 ± 11	96 ± 11	105 ± 6	105 ± 8	102 ± 8	96 ± 7	102 ± 7	94 ± 6	91 ± 8	92 ± 7
	4-Heptylphenol	6.35	8.47	5.39	7.22	6.98	6.61	4.68	9.35	8.91	9.30	5.25	10.34	5.84	9.37	6.80	7.44
		101 ± 9	91 ± 8	103 ± 9	99 ± 8	104 ± 7	93 ± 6	100 ± 7	105 ± 8	98 ± 5	105 ± 7	91 ± 8	101 ± 9	106 ± 6	98 ± 5	96 ± 7	95 ± 6
	Bisphenol Z	6.62	9.55	6.28	8.15	7.19	8.45	5.85	6.78	6.13	7.24	7.34	9.35	9.82	8.50	5.29	11.32
		98 ± 10	99 ± 9	99 ± 4	103 ± 10	105 ± 8	92 ± 5	91 ± 6	103 ± 6	106 ± 6	103 ± 9	100 ± 8	94 ± 7	98 ± 8	101 ± 10	91 ± 8	91 ± 8
	Pentachlorophenol	8.85	7.75	7.41	6.88	8.65	6.30	10.32	8.15	7.17	7.48	7.91	11.02	5.82	8.36	6.38	10.94
		105 ± 6	95 ± 8	91 ± 9	94 ± 7	105 ± 7	104 ± 7	105 ± 7	101 ± 6	106 ± 6	93 ± 6	100 ± 9	101 ± 10	98 ± 7	93 ± 6	99 ± 9	96 ± 8
	Triclosan	9.13	5.64	8.23	6.13	7.18	9.40	8.07	7.50	5.55	10.10	6.16	7.22	8.15	10.13	5.24	5.78
		94 ± 7	105 ± 7	105 ± 7	98 ± 10	97 ± 8	105 ± 6	96 ± 9	98 ± 6	104 ± 9	98 ± 8	103 ± 9	102 ± 6	104 ± 8	93 ± 7	105 ± 8	99 ± 8

^a RSD: relative standard deviation (n=12) for 500 and 2000 ng/kg (p< 0.05), interday. Percent recoveries (% ± SD, n = 3) of analytes spiked to milk and dairy product samples (500 – 2 000 ng/L)

Table S2. Results obtained in the matrix effect studies of the different types of dairy product samples

		Matrix Effect (%) ^a							
		Milk	Yoghourt	Cheese	Milkshake	Cream	Custard	Butter	Margarine
Parabens	Methylparaben	1.03 (3%)	0.98 (-2%)	0.94 (-6%)	1.06 (6%)	0.95 (-5%)	0.97 (-3%)	1.02 (2%)	0.92 (-7%)
	Ethylparaben	0.98 (-2%)	1.04 (4%)	1.05 (5%)	0.93 (-7%)	0.99 (-1%)	1.03 (3%)	0.95 (-5%)	0.93 (-7%)
	Isopropylparaben	1.07 (7%)	0.98 (-2%)	0.94 (-6%)	1.05 (5%)	0.95 (-5%)	1.03 (3%)	1.05 (5%)	0.93 (-7%)
	Butylparaben	1.02 (2%)	0.92 (-8%)	0.99 (-1%)	1.01 (1%)	0.95 (-5%)	0.90 (-10%)	1.01 (1%)	1.02 (2%)
	Propylparaben	1.02 (2%)	0.92 (-7%)	0.95 (-5%)	0.90 (-10%)	1.06 (6%)	0.92 (-8%)	0.90 (-10%)	0.98 (-2%)
	Isobutylparaben	0.95 (-5%)	0.93 (-7%)	0.93 (-7%)	1.04 (4%)	0.93 (-7%)	1.04 (4%)	1.04 (4%)	0.99 (-1%)
	Benzylparaben	0.99 (-1%)	1.05 (5%)	1.03 (3%)	0.98 (-2%)	1.03 (3%)	0.98 (-2%)	0.99 (-1%)	1.05 (5%)
Phenols	Pentylphenol	1.02 (2%)	0.95 (-5%)	0.90 (-10%)	0.98 (-2%)	1.02 (2%)	0.92 (-7%)	0.94 (-6%)	1.04 (4%)
	4-Chloro-3-methylphenol	1.07 (7%)	1.06 (6%)	0.92 (-8%)	0.93 (-7%)	0.93 (-7%)	0.97 (-3%)	0.90 (-10%)	1.04 (4%)
	4-tert-Octylphenol	0.95 (-5%)	0.93 (-7%)	0.95 (-5%)	0.93 (-7%)	1.02 (2%)	0.93 (-7%)	0.95 (-5%)	0.93 (-7%)
	Phenol	0.99 (-1%)	1.03 (3%)	0.95 (-5%)	1.04 (4%)	1.05 (5%)	1.06 (6%)	1.04 (4%)	1.02 (2%)
	4-tert-Butylphenol	0.95 (-5%)	1.03 (3%)	0.90 (-10%)	0.92 (-7%)	0.96 (-4%)	0.95 (-5%)	0.98 (-2%)	0.99 (-1%)
	Nonylphenol	0.98 (-2%)	0.93 (-7%)	0.97 (-3%)	0.93 (-7%)	0.94 (-6%)	1.03 (3%)	0.99 (-1%)	1.03 (3%)
	2-tert-Butylphenol	0.98 (-2%)	1.02 (2%)	0.93 (-7%)	1.02 (2%)	0.93 (-7%)	0.93 (-7%)	0.93 (-7%)	1.06 (6%)
	Bisphenol F	0.94 (-6%)	1.03 (3%)	0.99 (-1%)	1.03 (3%)	0.99 (-1%)	0.93 (-7%)	1.02 (2%)	1.10 (10%)
	3,4-Dimethylphenol	0.93 (-7%)	1.06 (6%)	1.02 (2%)	0.92 (-7%)	1.02 (2%)	0.93 (-7%)	0.92 (-8%)	0.92 (-8%)
	2,5-Dimethylphenol	0.98 (-2%)	0.92 (-7%)	0.93 (-7%)	1.04 (4%)	0.90 (-10%)	0.98 (-2%)	0.93 (-7%)	0.95 (-5%)
	4-Phenylphenol	0.95 (-5%)	1.04 (4%)	0.90 (-10%)	1.04 (4%)	1.02 (2%)	0.94 (-6%)	1.06 (6%)	0.95 (-5%)
	2-Phenylphenol	1.01 (1%)	1.04 (4%)	0.95 (-5%)	0.93 (-7%)	0.99 (-1%)	1.05 (5%)	1.05 (5%)	0.93 (-7%)
	Bisphenol S	0.99 (-1%)	0.93 (-7%)	1.05 (5%)	1.06 (6%)	0.98 (-2%)	0.92 (-7%)	1.06 (6%)	0.93 (-7%)
	4-Chlorophenol	1.04 (4%)	1.02 (2%)	0.96 (-4%)	0.95 (-5%)	0.99 (-1%)	1.05 (5%)	0.95 (-5%)	1.02 (2%)
	4-Hexylphenol	1.05 (5%)	1.01 (1%)	1.06 (6%)	0.96 (-6%)	0.92 (-8%)	1.02 (2%)	0.93 (-7%)	0.99 (-1%)
	Bisphenol A	1.02 (2%)	0.94 (-6%)	1.04 (4%)	1.01 (1%)	1.06 (6%)	0.98 (-2%)	1.08 (8%)	0.95 (-5%)
	Bisphenol B	0.99 (-1%)	1.05 (5%)	0.95 (-5%)	0.99 (-1%)	0.95 (-5%)	0.99 (-1%)	1.02 (2%)	0.95 (-5%)
	4-Heptylphenol	0.94 (-6%)	1.04 (4%)	1.05 (5%)	1.04 (4%)	0.98 (-2%)	1.02 (2%)	0.93 (-7%)	1.04 (4%)
	Bisphenol Z	0.92 (8%)	1.01 (1%)	1.08 (8%)	0.92 (-8%)	1.08 (8%)	1.05 (5%)	0.92 (-8%)	0.93 (-7%)
	Pentachlorophenol	0.90 (-10%)	1.02 (2%)	1.08 (8%)	1.02 (2%)	0.91 (-9%)	0.98 (-2%)	0.92 (-7%)	0.95 (-5%)
	Triclosan	1.02 (2%)	1.06 (6%)	0.91 (-9%)	1.06 (6%)	1.01 (1%)	0.93 (-7%)	1.06 (6%)	1.06 (6%)



Step	Valve 1	Valve 2	Valve 3
1. ERM-Lipid sorbent conditioning	L	I	I
2. Oasis PRiME HLB sorbent conditioning	I	L	I
3. Introduction of sample or standards	L	L	L
4. Oasis PRiME HLB drying	L	I	L
5. Elution of analytes	L	I	I

Figure S1. Continuous-flow system for the solid-phase with three injection valves and two serially columns.







Figure S2. Extract ion chromatograms for four analytes (benzylparaben, bisphenol A, triclosan and 4-heptylphenol). In one chromatogram, the analyte is spiked in the solvent, and in the other the analyte is spiked in a milk sample.