

Table S1. Primers used in the qPCR analyses in vitro.

Gene	Forward (5' to 3')	Reverse (5' to 3')
CD64	CTTCTCCTTCTATGTGGGCAGT	GCTACCTCGCACCCAGTATGAT
CD163	TCACTCCTGGGCTGCACGTAAAC	GATGTTATTTGCCATACAGGAGAATTG
IL-12	AAGATGAAGGAGACAGAG	CATTGGACTTCGGTAGAT
IL-10	GTGGAGCAGGTGAAGAGTGA	TTCATGGCCTTGTAGACACCT
GAPDH	ACAATGAATACGGCTACAG	GGTCCAGGGTTTCTTACT

Table S2. Primers used in the qPCR analyses in vivo.

Gene	Forward (5' to 3')	Reverse (5' to 3')
TNF- α	ACCTTGTTGCCTCCTCTT	GTTCACTGATGTAGCGACAG
IL-6	AAATGAGAAAAGAGTTGTG	TTTGTATCTCTGGAAGTTT
IL-1 β	GATACCACTCCCAACAGA	GCCATTGCACAACCTCTTT
F4/80	TGTCTGAAGATTCTCAAAACATGGA	TGGAGCTTCATAGTTGTAAGGCA
MMP-9	CGATTCCAAACCTTCAAA	GCAAGTCTTCAGAGTAGT
TIMP-1	AAGATGACTAAGATGCTAA	GATGAGAAACTCTTCACT
GAPDH	ACAATGAATACGGCTACAG	GGTCCAGGGTTTCTTACT