

Supplementary Materials




Figure S1. Purification of p67 Δ TM from HEK293T cells (a) and BLV Gag from *E. coli* (b). Samples taken before (cell media and cell lysate), during (flow-through, wash 1 and wash 2) and after (eluate and concentrated eluate) purification of His-tagged p67 Δ TM and BLV Gag through a cobalt-agarose column were used for SDS PAGE and Coomassie-blue staining. (c) Purified p67 Δ TM was diluted in TBS from 1:5 up to 1:500 and subjected to SDS PAGE and Coomassie-blue staining as undiluted purified sample results in overloading, as also seen in (a).