

Supplementary figure 1.



\* Correct sequence.



\* Correct sequence.

Supplementary Figure 1. Colony PCR to check proper cloning of pet32a-PirA<sup>vp</sup> and pet28b-PirB<sup>vp</sup>. Checking presence of inserts by amplifying the cloning site using gene-specific primers. Band size 336 bp and 1317 bp were observed, which correspond to PirA<sup>vp</sup> and PirB<sup>vp</sup>, respectively. One clone from each gene was sent for sequencing to further check the identity of the gene inserted in the plasmid.