



Supplementary Materials

Nanomaterial Characterization in Complex Media—Guidance and Application

Yves Uwe Hachenberger ^{1,†}, Daniel Rosenkranz ^{2,†}, Charlotte Kromer ¹, Benjamin Christoph Krause ^{1,†}, Nadine Dreiaek ¹, Fabian Lukas Kriegel ^{1,†}, Ekaterina Koz'menko ^{1,†}, Harald Jungnickel ¹, Jutta Tentschert ¹, Frank Stefan Bierkandt ^{1,*}, Peter Laux ^{1,*}, Ulrich Panne ³ and Andreas Luch ¹

¹ Department of Chemical & Product Safety, German Federal Institute for Risk Assessment (BfR), Max-Dohrn-Strasse 8-10, 10589 Berlin, Germany

² Institute for Clinical Chemistry and Laboratory Medicine, Klinikum Oldenburg AöR, Rahel-Straus-Straße 10, 26133 Oldenburg, Germany

³ Federal Institute for Materials Research and Testing (BAM), Richard-Willstätter-Strasse 11, 12489 Berlin, Germany

* Correspondence: frank.bierkandt@bfr.bund.de (F.S.B.); peter.laux@bfr.bund.de (P.L.)

† Former Employee of BfR or/and BAM.




Figure S1. Time-resolved DI spICP-MS analysis of CuO NP dispersions (1 µg/L) in aqueous media.

Time-resolved DI spICP-MS Measurements of CuO NP in BSA




Figure S2. Time-resolved DI spICP-MS analysis of CuO NP dispersions (1 $\mu\text{g/L}$) in BSA.

Time-resolved DI spICP-MS Measurements of CuO NP in NaHCO_3




Figure S3. Time-resolved DI spICP-MS analysis of CuO NP dispersions (1 $\mu\text{g/L}$) in NaHCO_3 .

Time-resolved DI spICP-MS Measurements of CuO NP in ALF




Figure S4. Time-resolved DI spICP-MS analysis of CuO NP dispersions (1 $\mu\text{g/L}$) in ALF.

Time-resolved DI spICP-MS Measurements of CuO NP in GS




Figure S5. Time-resolved DI spICP-MS analysis of CuO NP dispersions (1 $\mu\text{g/L}$) in GS.

Time-resolved DI spICP-MS Measurements of CuO NP in EGS




Figure S6. Time-resolved DI spICP-MS analysis of CuO NP dispersions (1 $\mu\text{g/L}$) in EGS.

Time-resolved DI spICP-MS Measurements of CuO NP in RPMI




Figure S7. Time-resolved DI spICP-MS analysis of CuO NP dispersions (1 $\mu\text{g/L}$) in RPMI.